

# Package ‘OPI’

January 24, 2024

**Type** Package

**Title** Open Perimetry Interface

**Version** 2.11.2

**Date** 2024-01-24

**Author** Andrew Turpin [cre, aut, cph],  
David Lawson [ctb, cph],  
Matthias Muller [ctb],  
Jonathan Dennis [ctb, cph],  
Astrid Zeman [ctb],  
Ivan Marin-Franch [ctb]

**Maintainer** Andrew Turpin <andrew.turpin@lei.org.au>

**Description** Implementation of the Open Perimetry Interface (OPI) for simulating and controlling visual field machines using R. The OPI is a standard for interfacing with visual field testing machines (perimeters) first started as an open source project with support of Haag-Streit in 2010. It specifies basic functions that allow many visual field tests to be constructed. As of February 2022 it is fully implemented on the Haag-Streit Octopus 900 with partial implementations on the Centervue Compass, Kowa AP 7000, Android phones and the CrewT IMO. It also has a cousin: the R package 'visual-Fields', which has tools for analysing and manipulating visual field data.

**License** Apache License (== 2.0)

**URL** <https://people.eng.unimelb.edu.au/aturpin/opi/index.html>

**Depends** methods, Rfast, abind, grDevices, openssl

**LazyData** true

**Encoding** UTF-8

**RoxygenNote** 7.3.1

**Collate** OPI-package.r opi.r compassClient.r displayClient.r  
daydreamClient.r dbToCdr.r fourTwo.r full\_threshold.r  
imoClient.r kowaAP7000Client.r mocs.r octopus600.r  
octopus900Client.r pix2deg.r phoneHMD.r QUESTP.r simDisplay.r  
simG.r simH.r simH\_RT.r simNo.r simYes.r zest.r  
data-RtDbUnits.r data-RtSigmaUnits.r opiKineticStimulus.r  
opiStaticStimulus.r opiTemporalStimulus.r

**NeedsCompilation** no

**Repository** CRAN

**Date/Publication** 2024-01-24 09:12:46 UTC

## R topics documented:

.OpiEnv . . . . .	2
cdTodb . . . . .	3
chooseOpi . . . . .	3
dbTocd . . . . .	5
degTopix . . . . .	5
fourTwo.start . . . . .	6
FT . . . . .	8
MOCS . . . . .	11
opi.implementations . . . . .	16
opiClose . . . . .	16
opiGetParams . . . . .	18
opiInitialize . . . . .	19
opiKineticStimulus . . . . .	26
opiPresent . . . . .	28
opiQueryDevice . . . . .	35
opiSetBackground . . . . .	37
opiStaticStimulus . . . . .	44
opiTemporalStimulus . . . . .	46
pixTodeg . . . . .	47
QUESTP . . . . .	48
RtDbUnits . . . . .	56
RtSigmaUnits . . . . .	57
ZEST . . . . .	58
<b>Index</b>	<b>63</b>

---

.OpiEnv

*Environment holding the state of the OPI*

---

### Description

Holds the chosen machine and any parameters forming the state of the machine.

### Usage

.OpiEnv

### Format

An object of class environment of length 13.

---

cdTodb	<i>Convert cd/m<sup>2</sup> to perimetric dB.</i>
--------	---

---

**Description**

Given a value in  $\text{cd/m}^2$ , return the equivalent dB value. Default is to use HFA units, so maximum stimulus is 10000 apostilbs.

**Usage**

```
cdTodb(cd, maxStim = 10000/pi)
```

**Arguments**

cd	Value to convert to dB in $\text{cd/m}^2$ . Must be $> 0$ .
maxStim	Stimulus value for 0dB in $\text{cd/m}^2$ .

**Value**

A dB value for cd  $\text{cd/m}^2$ .

**Examples**

```
# candela to decibels
dB <- cdTodb(10000/pi) # 0 dB
dB <- cdTodb(1000/pi) # 10 dB
dB <- cdTodb(100/pi) # 20 dB
dB <- cdTodb(10/pi) # 30 dB
dB <- cdTodb(1/pi) # 40 dB
dB <- cdTodb(0.1/pi) # 50 dB
```

---

chooseOpi	<i>Choose an implementation of the OPI</i>
-----------	--

---

**Description**

Chooses an implementation of the OPI to use

**Usage**

```
chooseOpi(opiImplementation)

chooseOPI(opiImplementation)
```

**Arguments**

opiImplementation

A character string that is one of the following.

- "SimNo" for a simulator that always doesn't see.
- "SimYes" for a simulator that always does see.
- "SimHenson" for a simulator that uses a cumulative gaussian psychometric function with standard deviation according to Henson et al (2000) where variability increases as true threshold (Humphrey dB) value decreases.
- "SimHensonRT" as for SimHenson, but response times in ms are sampled from a supplied response time data set for each true positive response.
- "SimGaussian" for a simulator that uses a cumulative gaussian psychometric function with standard deviation supplied in opiInitialize().
- "Octopus900" for interfacing with the Octopus 900.
- "Octopus900F310" for interfacing with the Octopus 900 using Logitech F310 controller.
- "Octopus600" for interfacing with the Octopus 600.
- "HEP" not working so well in HEPs.
- "KowaAP7000" for interfacing with Kowa AP-7000.
- "Imo" for interfacing with CrewT's Imo head mounted perimeter.
- "DayDream" for interfacing with an Android phone in a Google Daydream.
- "Display" for interfacing with a shiny plot area on the current machine.
- "PhoneHMD" for interfacing with phones using VR. At the moment, only Android compatible phones are working. The VR headset must be compatible with Cardboard
- NULL print a list of available OPI implementations.

**Value**

Returns TRUE if successful, FALSE otherwise.

**References**

David B. Henson, Shaila Chaudry, Paul H. Artes, E. Brian Faragher, and Alec Ansons. Response Variability in the Visual Field: Comparison of Optic Neuritis, Glaucoma, Ocular Hypertension, and Normal Eyes. *Investigative Ophthalmology & Visual Science*, February 2000, Vol. 41(2).

A.M. McKendrick, J. Denniss and A. Turpin. "Response times across the visual field: empirical observations and application to threshold determination". *Vision Research*, 101, 2014.

A. Turpin, P.H. Artes and A.M. McKendrick. "The Open Perimetry Interface: An enabling tool for clinical visual psychophysics", *Journal of Vision* 12(11) 2012.

**Examples**

```
if(!chooseOpi("SimHenson"))
  warnings()
```

---

dbToCd	<i>Convert perimetric dB to cd/m<sup>2</sup></i>
--------	--

---

**Description**

Given a value in dB, return the cd/m<sup>2</sup> equivalent. Default is to use HFA units, so maximum stimulus is 10000 apostilbs.

**Usage**

```
dbToCd(db, maxStim = 10000/pi)
```

**Arguments**

db	Value to convert to cd/m <sup>2</sup> .
maxStim	Stimulus value for 0dB in cd/m <sup>2</sup> .

**Value**

cd/m<sup>2</sup> value for db dB.

**Examples**

```
# decibels to candela
cd <- dbToCd(0) # 10000/pi
cd <- dbToCd(10) # 1000/pi
cd <- dbToCd(20) # 100/pi
cd <- dbToCd(30) # 10/pi
cd <- dbToCd(40) # 1/pi
```

---

degToPix	<i>Convert degrees to pixels for machine 'machine'</i>
----------	--

---

**Description**

Convert degrees to pixels for machine 'machine'

**Usage**

```
degToPix(xy, machine = "compass")
```

**Arguments**

xy	a 2 element vector c(x,y) where x and y are in pixels
machine	"compass" or ...?

**Value**

xy converted to pixels (top-left is (0,0)) for the machine or NA if machine is unknown

**Examples**

```
degTopix(c(0, 0), machine="compass") # c(960, 960) pixels
degTopix(c(-15, 2)) # c(495, 898) pixels
```

---

fourTwo.start	<i>4-2 Staircase</i>
---------------	----------------------

---

**Description**

fourTwo is a 4-2 dB staircase beginning at level est terminating after two reversals. The final estimate is the average of the last two presentations. It also terminates if the minStimulus is not seen twice, or the maxStimulus is seen twice.

**Usage**

```
fourTwo.start(est = 25, instRange = c(0, 40), verbose = FALSE, makeStim, ...)
fourTwo.step(state, nextStim = NULL)
fourTwo.stop(state)
fourTwo.final(state)
```

**Arguments**

est	Starting estimate in dB
instRange	Dynamic range of the instrument c(min,max) in dB
verbose	True if you want each presentation printed
makeStim	A function that takes a dB value and numPresentations and returns an OPI datatype ready for passing to opiPresent
...	Extra parameters to pass to the opiPresent function
state	Current state of the fourTwo returned by fourTwo.start and fourTwo.step
nextStim	A valid object for opiPresent to use as its nextStim.

**Details**

This is an implementation of a 4-2 1-up 1-down staircase. The initial staircase starts at est and proceeds in steps of 4 dB until the first reversal, and 2dB until the next reversal. The mean of the last two presentations is taken as the threshold value. Note this function will repeatedly call opiPresent for a stimulus until opiPresent returns NULL (ie no error occurred). If more than one fourTwo is to be interleaved (for example, testing multiple locations), then the fourTwo.start, fourTwo.step, fourTwo.stop and fourTwo.final calls can maintain the state of the fourTwo after each presentation, and should be used. See examples below.

**Value****Multilple locations:**

fourTwo.start returns a list that can be passed to fourTwo.step, fourTwo.stop, and fourTwo.final.

It represents the state of a fourTwo at a single location at a point in time and contains the following.

\* name, fourTwo. \* startingEstimate=est, input param. \* currentLevel, the next stimulus to present. \* minStimulus=instRange[1], input param. \* maxStimulus=instRange[2], input param. \* makeStim, input param. \* lastSeen, the last seen stimulus. \* lastResponse, the last response given. \* stairResult, The final result if finished (initially NA). \* finished, "Not" if staircase has not finished, or one of "Rev" (finished due to 2 reversals), "Max" (finished due to 2 maxStimulus seen), "Min" (finished due to 2 minStimulus not seen). \* verbose, number of reversals so far. \* numberOfReversals, number of reversals so far. \* currSeenLimit, number of times maxStimulus has been seen. \* currNotSeenLimit, number of times minStimulus not seen. \* numPresentations, number of presentations so far. \* stimuli, vector of stimuli shown at each call to fourTwo.step. \* responses, vector of responses received (1 seen, 0 not) received at each call to fourTwo.step. \* responseTimes, vector of response times received at each call to fourTwo.step. \* opiParams=list(...), input param

fourTwo.step returns a list containing \* state, the new state after presenting a stimuli and getting a response. \* resp, the return from the opiPresent call that was made.

fourTwo.stop returns TRUE if the staircase is finished (2 reversals, or maxStimulus is seen twice or minStimulus is not seen twice).

fourTwo.final returns the final estimate of threshold (mean of last two reversals). This issues a warning if called before the staircase has finished, but still returns a value.

**See Also**

[dbTocd](#), [opiPresent](#), [FT](#)

**Examples**

```
# Stimulus is Size III white-on-white as in the HFA
makeStim <- function(db, n) {
  s <- list(x=9, y=9, level=dbTocd(db), size=0.43, color="white",
           duration=200, responseWindow=1500)
  class(s) <- "opiStaticStimulus"
  return(s)
}
chooseOpi("SimHenson")
if (!is.null(opiInitialize(type="C", cap=6)))
  stop("opiInitialize failed")

#####
# This section is for multiple fourTwos
#####
makeStimHelper <- function(db,n, x, y) { # returns a function of (db,n)
  ff <- function(db, n) db+n
  body(ff) <- substitute({
    s <- list(x=x, y=y, level=dbTocd(db), size=0.43, color="white",
             duration=200, responseWindow=1500)
    class(s) <- "opiStaticStimulus"
    return(s)}, list(x=x,y=y))
}
```

```

    return(ff)
  }
  # List of (x, y, true threshold) triples
  locations <- list(c(9,9,30), c(-9,-9,32), c(9,-9,31), c(-9,9,33))

  # Setup starting states for each location
  states <- lapply(locations, function(loc) {
    fourTwo.start(makeStim=makeStimHelper(db,n,loc[1],loc[2]),
                  tt=loc[3], fpr=0.03, fnr=0.01)})

  # Loop through until all states are "stop"
  while(!all(st <- unlist(lapply(states, fourTwo.stop)))) {
    i <- which(!st) # choose a random,
    i <- i[runif(1, min=1, max=length(i))] # unstopped state
    r <- fourTwo.step(states[[i]]) # step it
    states[[i]] <- r$state # update the states
  }

  finals <- lapply(states, fourTwo.final) # get final estimates of threshold
  for(i in 1:length(locations)) {
    cat(sprintf("Location (%+2d,%+2d) ",locations[[i]][1], locations[[i]][2]))
    cat(sprintf("has threshold %4.2f\n", finals[[i]]))
  }

  if (!is.null(opiClose()))
    warning("opiClose() failed")

```

---

 FT

*Full Threshold*


---

## Description

FT begins with a 4-2dB staircase beginning at level est. If the final estimate (last seen) is more than 4dB away from est, a second 4-2 staircase is completed beginning at the estimate returned from the first

## Usage

```
FT(est = 25, instRange = c(0, 40), verbose = FALSE, makeStim, ...)
```

```
FT.start(est = 25, instRange = c(0, 40), makeStim, ...)
```

```
FT.step(state, nextStim = NULL)
```

```
FT.stop(state)
```

```
FT.final(state)
```



**Arguments**

<code>est</code>	Starting estimate in dB
<code>instRange</code>	Dynamic range of the instrument c(min,max) in dB
<code>verbose</code>	True if you want each presentation printed
<code>makeStim</code>	A function that takes a dB value and numPresentations and returns an OPI datatype ready for passing to opiPresent
<code>...</code>	Extra parameters to pass to the opiPresent function
<code>state</code>	Current state of the FT returned by FT.start and FT.step
<code>nextStim</code>	A valid object for opiPresent to use as its nextStim

**Details**

This is an implementation of a 4-2 1-up 1-down staircase as implemented in the first Humphrey Field Analyzer. The initial staircase starts at `est` and proceeds in steps of 4 dB until the first reversal, and 2dB until the next reversal. The last seen stimulus is taken as the threshold value. If, after the first staircase, the threshold is more than 4 dB away from the starting point, then a second staircase is initiated with a starting point equal to the threshold found with the first staircase.

Note this function will repeatedly call `opiPresent` for a stimulus until `opiPresent` returns NULL (ie no error occurred)

If more than one FT is to be interleaved (for example, testing multiple locations), then the `FT.start`, `FT.step`, `FT.stop` and `FT.final` calls can maintain the state of the FT after each presentation, and should be used. If only a single FT is required, then the simpler FT can be used. See examples below

**Value****Single location:**

Returns a list containing `* npres`, total number of presentations. `* respSeq`, response sequence stored as a list of (seen,dB) pairs. `* first`, first staircase estimate in dB. `* final`, final threshold estimate in dB.

**Multilple locations:**

`FT.start` returns a list that can be passed to `FT.step`, `FT.stop`, and `FT.final`. It represents the state of a FT at a single location at a point in time and contains the following. `* name`, FT. `* startingEstimate=est`, input param. `* currentLevel`, the next stimulus to present. `* minStimulus=instRange[1]`, input param. `* maxStimulus=instRange[2]`, input param. `* makeStim`, input param. `* lastSeen`, the last seen stimulus. `* lastResponse`, the last response given. `* stairResult`, The final result if finished (initially NA). `* finished`, "Not" if staircase has not finished, or one of "Rev" (finished due to 2 reversals), "Max" (finished due to 2 maxStimulus seen), "Min" (finished due to 2 minStimulus not seen). `* verbose`, number of reversals so far. `* numberOfReversals`, number of reversals so far. `* currSeenLimit`, number of times maxStimulus has been seen. `* currNotSeenLimit`, number of times minStimulus not seen. `* numPresentations`, number of presentations so far. `* stimuli`, vector of stimuli shown at each call to `FT.step`. `* responses`, vector of responses received (1 seen, 0 not) received at each call to `FT.step`. `* responseTimes`, vector of response times received at each call to `FT.step`. `*`

`opiParams=list(...)`, input param \* finished, TRUE if staircase has finished (2 reversals, or max/min seen/not-seen twice).

`FT.step` returns a list containing

- `state`, the new state after presenting a stimuli and getting a response.
- `resp`, the return from the `opiPresent` call that was made.

`FT.stop` returns TRUE if the first staircase has had 2 reversals, or `maxStimulus` is seen twice or `minStimulus` is not seen twice and the final estimate is within 4 dB of the starting stimulus. Returns TRUE if the second staircase has had 2 reversals, or `maxStimulus` is seen twice or `minStimulus` is not seen twice

`FT.final` returns the final estimate of threshold based on state, which is the last seen in the second staircase, if it ran, or the first staircase otherwise

`FT.final.details` returns a list containing

- `final`, the final threshold.
- `first`, the threshold determined by the first staircase (might be different from final).
- `stopReason`, either Reversals, Max, or Min which are the three ways in which FT can terminate.
- `np`, number of presentation for the whole procedure (including both staircases if run).

## References

A. Turpin, P.H. Artes and A.M. McKendrick. "The Open Perimetry Interface: An enabling tool for clinical visual psychophysics", *Journal of Vision* 12(11) 2012.

H. Bebie, F. Fankhauser and J. Spahr. "Static perimetry: strategies", *Acta Ophthalmology* 54 1976.

C.A. Johnson, B.C. Chauhan, and L.R. Shapiro. "Properties of staircase procedures for estimating thresholds in automated perimetry", *Investagative Ophthalmology and Vision Science* 33 1993.

## See Also

[dbToCd](#), [opiPresent](#), [fourTwo.start](#)

## Examples

```
# Stimulus is Size III white-on-white as in the HFA
makeStim <- function(db, n) {
  s <- list(x=9, y=9, level=dbToCd(db), size=0.43, color="white",
           duration=200, responseWindow=1500)
  class(s) <- "opiStaticStimulus"
  return(s)
}
chooseOpi("SimHenson")
if (!is.null(opiInitialize(type="C", cap=6)))
  stop("opiInitialize failed")

result <- FT(makeStim=makeStim, tt=30, fpr=0.15, fnr=0.01)
if (!is.null(opiClose()))
  warning("opiClose() failed")
```

```
#####
```

```

# This section is for multiple FTs
#####
makeStimHelper <- function(db,n, x, y) { # returns a function of (db,n)
  ff <- function(db, n) db+n
  body(ff) <- substitute({
    s <- list(x=x, y=y, level=dbTocd(db), size=0.43, color="white",
             duration=200, responseWindow=1500)
    class(s) <- "opiStaticStimulus"
    return(s)
  }, list(x=x,y=y))
  return(ff)
}

# List of (x, y, true threshold) triples
locations <- list(c(9,9,30), c(-9,-9,32), c(9,-9,31), c(-9,9,33))
# Setup starting states for each location
states <- lapply(locations, function(loc) {
  FT.start(makeStim=makeStimHelper(db,n,loc[1],loc[2]),
           tt=loc[3], fpr=0.03, fnr=0.01)})

# Loop through until all states are "stop"
while(!all(st <- unlist(lapply(states, FT.stop)))) {
  i <- which(!st) # choose a random,
  i <- i[runif(1, min=1, max=length(i))] # unstopped state
  r <- FT.step(states[[i]]) # step it
  states[[i]] <- r$state # update the states
}

finals <- lapply(states, FT.final) # get final estimates of threshold
for(i in 1:length(locations)) {
  cat(sprintf("Location (%+2d,%+2d) ",locations[[i]][1], locations[[i]][2]))
  cat(sprintf("has threshold %4.2f\n", finals[[i]]))
}

if(!is.null(opiClose()))
  warning("opiClose() failed")

```

---

MOCS

*Method of Constant Stimuli (MOCS)*


---

### Description

MOCS performs either a yes/no or n-interval-forced-choice Method of Constant Stimuli test

### Usage

```

MOCS(
  params = NA,
  order = "random",
  responseWindowMeth = "constant",

```

```

responseFloor = 1500,
responseHistory = 5,
keyHandler = function(correct, ret) return(list(seen = TRUE, time = 0, err = NULL)),
interStimMin = 200,
interStimMax = 500,
beep_function,
makeStim,
stim_print,
...
)

```

### Arguments

params	<p>A matrix where each row is <math>x\ y\ i\ n\ correct\_n\ l11\ l12\ \dots\ l1m</math> where</p> <ul style="list-style-type: none"> <li>• <math>x</math> is X coordinate of location.</li> <li>• <math>y</math> is Y coordinate of location.</li> <li>• <math>i</math> is a location number (assigned by caller).</li> <li>• <math>n</math> is Number of times this location/luminance(s) should be repeated.</li> <li>• <math>correct\_n</math> is the index <math>i</math> of the luminance level (<math>l1i</math>) that should be treated as a “correct” response (the correct interval). For a standard MOCS, this will be 1; for a 2AFC, this will be 1 or 2. This number will be in the range <math>[1, m]</math>.</li> <li>• <math>l1i</math> is the <math>i</math>'th luminance level to be used at this location for interval <math>i</math> of the presentation in <math>cd/m^2</math>. For a standard MOCS, <math>i=1</math>, and the params matrix will have 5 columns. For a 2AFC, there will be two <math>l1i</math>'s, and params will have 6 columns.</li> </ul>
order	<p>Control the order in which the stimuli are presented.</p> <ul style="list-style-type: none"> <li>• "random" Randomise the order of trials/locations.</li> <li>• "fixed" Present each row of params in order of <math>1:nrow(params)</math>, ignoring the <math>n</math> (4th) column in params.</li> </ul>
responseWindowMeth	<p>Control time perimeter waits for response.</p> <ul style="list-style-type: none"> <li>• "speed" After an average of the last speedHistory response times, with a minimum of responseFloor. Initially responseFloor.</li> <li>• "constant" Always use responseFloor.</li> <li>• "forceKey" Wait for a keyboard input.</li> </ul>
responseFloor	Minimum response window (for any responseWindowMeth except "forceKey").
responseHistory	Number of past yeses to average to get response window (only used if responseWindowMeth is "speed").
keyHandler	Function to get a keyboard input and returns as for <code>opiPresent</code> : <code>list(err={NULL msg}, seen={TRUE FALSE}, time)</code> . The parameters passed to the function are the correct interval number (column 4 of params), and the result of <code>opiPresent</code> . See Examples.
interStimMin	Regardless of response, wait <code>runif(interStimMin, interStimMax)</code> ms.

<code>interStimMax</code>	Regardless of response, wait <code>runif(interStimMin, interStimMax)</code> ms.
<code>beep_function</code>	A function that takes the string 'correct', the string 'incorrect', or a stimulus number and plays an appropriate sound. See examples.
<code>makeStim</code>	A helper function to take a row of params and a response window length in ms, and create a list of OPI stimuli types for passing to <code>opiPresent</code> . This may include a <code>checkFixationOK</code> function. See Example.
<code>stim_print</code>	A function that takes an <code>opiStaticStimulus</code> and return list from <code>opiPresent</code> and returns a string to print for each presentation. It is called immediately after each <code>opiPresent</code> , and the string is prepended with the (x,y) coordinates of the presentation and ends with a newline.
<code>...</code>	Extra parameters to pass to the <code>opiPresent</code> function.

### Details

Whether the test is yes/no or forced-choice is determined by the number of columns in params. The code simply presents all columns from 5 onwards and collects a response at the end. So if there is only 5 columns, it is a yes/no task. If there are 6 columns it is a 2-interval-forced-choice. Generally, an nIFC experiment has 4+n columns in params.

Note that when the order is "random", the number of trials in the test will be the sum of the 3rd column of params. When the order is "fixed", there is only one presentation per row, regardless of the value in the 3rd column of params.

If a response is received before the final trial in a nIFC experiment, it is ignored.

If the `checkFixationOK` function is present in a stimulus, then it is called after each presentation, and the result is "anded" with each stimulus in a trial to get a TRUE/FALSE for fixating on all stimuli in a trial.

### Value

Returns a data.frame with one row per stimulus copied from params with extra columns appended: `checkFixation` checks, and the return values from `opiPresent()` (see example). These last values will differ depending on which machine/simulation you are running (as chosen with `chooseOpi()`).

- column 1: x
- column 2: y
- column 3: location number
- column 4: number of times to repeat this stim
- column 5: correct stimulus index
- column 6: TRUE/FALSE was fixating for all presentations in this trial according to `checkFixationOK`
- column 7...: columns from params
- ...: columns from `opiPresent` return

### References

A. Turpin, P.H. Artes and A.M. McKendrick. "The Open Perimetry Interface: An enabling tool for clinical visual psychophysics", *Journal of Vision* 12(11) 2012.

**See Also**

[dbTocd](#), [opiPresent](#)

**Examples**

```
# For the Octopus 900
# Check if pupil centre is within 10 pixels of (160,140)
checkFixationOK <- function(ret) return(sqrt((ret$pupilX - 160)^2 + (ret$pupilY - 140)^2) < 10)

# Return a list of opi stim objects (list of class opiStaticStimulus) for each level (dB) in
# p[5:length(p)]. Each stim has responseWindow BETWEEN_FLASH_TIME, except the last which has
# rwin. This one assumes p is on old Octopus 900 dB scale (0dB == 4000 cd/m^2).
makeStim <- function(p, rwin) {
  BETWEEN_FLASH_TIME <- 750 # ms
  res <- NULL
  for(i in 5:length(p)) {
    s <- list(x=p[1], y=p[2], level=dbTocd(p[i],4000/pi), size=0.43, duration=200,
             responseWindow=ifelse(i < length(p), BETWEEN_FLASH_TIME, rwin),
             checkFixationOK=NULL)
    class(s) <- "opiStaticStimulus"
    res <- c(res, list(s))
  }
  return(res)
}

#####
# Read in a key press 'z' is correct==1, 'm' otherwise
# correct is either 1 or 2, whichever is the correct interval
#
# Return list(seen={TRUE|FALSE}, time=time, err=NULL)
# seen is TRUE if correct key pressed
#####
## Not run:
if (length(dir(".", "getKeyPress.py")) < 1)
  stop('Python script getKeyPress.py missing?')

## End(Not run)

keyHandler <- function(correct, ret) {
  return(list(seen=TRUE, time=0, err=NULL))
  ONE <- "b'z'"
  TWO <- "b'm'"
  time <- Sys.time()
  key <- 'q'
  while (key != ONE && key != TWO) {
    a <- system('python getKeyPress.py', intern=TRUE)
    key <- a # substr(a, nchar(a), nchar(a))
    print(paste('Key pressed: ',key,'from',a))
    if (key == "b'8'")
      stop('Key 8 pressed')
  }
  time <- Sys.time() - time
}
```

```

    if ((key == ONE && correct == 1) || (key == TWO && correct == 2))
      return(list(seen=TRUE, time=time, err=NULL))
    else
      return(list(seen=FALSE, time=time, err=NULL))
  }

#####
# Read in return value from opipresent with F310 controller.
# First param is correct, next is 1 for left button, 2 for right button
# Left button (LB) is correct for interval 1, RB for interval 2
# correct is either 1 or 2, whichever is the correct interval
#
# Return list(seen={TRUE|FALSE}, time=time, err=NULL)
# seen is TRUE if correct key pressed
#####
F310Handler <- function(correct, opiResult) {
  z <- opiResult$seen == correct
  opiResult$seen <- z
  return(opiResult)
}

#####
# 2 example beep_function
#####
## Not run:
require(beepr)
myBeep <- function(type='None') {
  if (type == 'correct') {
    beepr::beep(2) # coin noise
    Sys.sleep(0.5)
  }
  if (type == 'incorrect') {
    beepr::beep(1) # system("rundll32 user32.dll,MessageBeep -1") # system beep
    #Sys.sleep(0.0)
  }
}
require(audio)
myBeep <- function(type="None") {
  if (type == 'correct') {
    wait(audio::play(sin(1:10000/10)))
  }
  if (type == 'incorrect') {
    wait(audio::play(sin(1:10000/20)))
  }
}

## End(Not run)

#####
# An example stim_print function
#####
## Not run:
stim_print <- function(s, ret) {

```

```

        sprintf("%4.1f %2.0f", cdTodb(s$level, 10000/pi), ret$seen)
    }

    ## End(Not run)

```

---

opi.implementations     *FOR INTERNAL USE ONLY*

---

### Description

The method `opiDistributor` searches for the specific method of a general OPI operation, which depends on the OPI implementation selected with `chooseOpi`. It returns an error if no OPI implementation has been selected yet. A catalog of all specific methods are listed in `opi.implementations`.

### Usage

```

opi.implementations

opiDistributor(operation, ...)

```

### Arguments

<code>operation</code>	A general OPI operation of the following methods to: <code>opiInitialize</code> , <code>opiPresent</code> , <code>opiClose</code> , <code>opiSetBackground</code> , <code>opiQueryDevice</code>
<code>...</code>	other parameters to pass to the methods

### Format

`opi.implementations` is a list containing a catalog of all specific methods that are dependent on the OPI implementation selected with `chooseOpi`

---

`opiClose`     *Close using OPI*

---

### Description

Generic function for closing the chosen OPI implementation that is set with `chooseOpi()`



**Usage**

```
opiClose(...)  
  
compass.opiClose()  
  
display.opiClose()  
  
daydream.opiClose()  
  
imo.opiClose()  
  
kowaAP7000.opiClose()  
  
octo600.opiClose()  
  
octo900.opiClose()  
  
PhoneHMD.opiClose()  
  
simG.opiClose()
```

**Arguments**

...                   Implementation specific parameters. See details.

**Value**

Returns NULL if close succeeded, otherwise an implementation-dependent error.

**Compass:** Returns a list of `err`, which is an error code, and `fixations`, which is a matrix with three columns: `time` (same as `time_hw` in `opiPresent`), `x` (degrees relative to the centre of the image returned by `opiInitialise` - not the PRL), `y` (as for `x`), and one row per fixation.

**Display:** Shuts the display.

**Daydream:** DETAILS

**imo:** DETAILS

**KowaAP7000:** DETAILS

**Octopus600:**  
Always returns NULL

**Octopus900:**  
Returns NULL.

**PhoneHMD:**  
Closes the socket connection with the PhoneHMD

**SimGaussian:**  
DETAILS

**See Also**[chooseOpi](#)**Examples**

```
chooseOpi("SimGaussian")
if (!is.null(opiInitialize(sd=2)))
  stop("opiInitialize failed")
if (!is.null(opiClose()))
  stop("opiClose failed, which is very surprising!")
```

---

`opiGetParams`*Get OPI method parameters*

---

**Description**

Get parameters of OPI functions which depends on the implementation set with `chooseOPI()`

**Usage**

```
opiGetParams(method, ...)
```

**Arguments**

<code>method</code>	Method for which to get parameters and defaults.
<code>...</code>	Implementation specific parameters. See details.

**Value**

Returns a list of parameters and their default vlues of the method `method` depending on the OPI implementation selected with `chooseOPI()`.

**Examples**

```
chooseOpi("SimHenson")
opiGetParams("opiInitialize")
opiGetParams("opiPresent")
```

---

opiInitialize	<i>Initialize OPI</i>
---------------	-----------------------

---

**Description**

Generic function for initialization of the chosen OPI implementation that is set with chooseOpi()

**Usage**

```
opiInitialize(...)

opiInitialise(...)

compass.opiInitialize(ip = "192.168.1.2", port = 44965)

display.opiInitialize(
  width,
  height,
  ppi,
  viewdist,
  lut = seq(0, 400, length.out = 256)
)

daydream.opiInitialize(
  ip = "127.0.0.1",
  port = 50008,
  lut = seq(0, 400, length.out = 256),
  fovy = 90
)

imo.opiInitialize(
  ip = "localhost",
  port = 1234,
  ppi = 16,
  tracking = FALSE,
  tracktol = 2
)

kowaAP7000.opiInitialize(ip = "192.168.1.2", port = 44965)

octo600.opiInitialize(
  ipAddress = "",
  eye = "",
  pupilTracking = FALSE,
  pulsar = FALSE,
  eyeControl = 0
)
```

```

octo900.opiInitialize(
  serverPort = 50001,
  eyeSuiteSettingsLocation = "C:/ProgramData/Haag-Streit/EyeSuite/",
  eye = "",
  gazeFeed = "",
  bigWheel = FALSE,
  pres_buzzer = 0,
  resp_buzzer = 0,
  zero_dB_is_10000_asb = TRUE
)

PhoneHMD.opiInitialize(ip, port = 50008, lut = seq(0, 400, length.out = 256))

simG.opiInitialize(sd = 2, display = NA, maxStim = 10000/pi)

simH.opiInitialize(
  type = "C",
  A = -0.081,
  B = 3.27,
  cap = 6,
  display = NA,
  maxStim = 10000/pi
)

simH_RT.opiInitialize(
  type = "C",
  cap = 6,
  A = -0.081,
  B = 3.27,
  display = NA,
  maxStim = 10000/pi,
  rtData = NULL,
  rtFP = 1:1600
)

```

### Arguments

...	Implementation specific parameters. See details.
ip	IP address on which server is listening for PhoneHMD
port	Port number on which server is listening for PhoneHMD. Default is 50008
width	Width of the screen in pixels
height	Height of the screen in pixels
ppi	Pixels per inch of the display
viewdist	Viewing distance in cm
lut	Look up table mapping pixel values to cd/m2

fovy	Field of view in degrees in the y-axis. It is different depending on the device. For Daydream view, it is 90 degrees, for, Daydream view 2 is 100 degrees. Default is 90.
ppd	pixels size as in pixels per degree
tracking	tracking on or off
tracktol	tolerance during tracking in degrees of visual angle
ipAddress	IP address of Octopus 600 machine
eye	eye; "right" or "left" for "Octopus900", "Octopus600"
pupilTracking	pupil tracking
pulsar	DETAILS
eyeControl	DETAILS
serverPort	port number on which server is listening for "Octopus900"
eyeSuiteSettingsLocation	dir name containing EyeSuite settings for "Octopus900"
gazeFeed	NA or a folder name for "Octopus900"
bigWheel	FALSE (standard machine), TRUE for modified aperture wheel for "Octopus900"
pres_buzzer	0 (no buzzer), 1, 2, 3 (max volume) for "Octopus900"
resp_buzzer	0 (no buzzer), 1, 2, 3 (max volume) for "Octopus900"
zero_dB_is_10000_asb	Is 0 dB 10000 apostilb (TRUE) or 4000 (FALSE) for "Octopus900"
sd	standard deviation for the Gaussian
display	Dimensions of plot area (-x,+x,-y,+y) to display stim. No display if NULL. For "SimHenson", "SimHensonRT", "SimGaussian", "SimNo", "SimYes"
maxStim	Maximum stimulus value in cd/m <sup>2</sup> used for db <-> cd/m <sup>2</sup> conversions for "SimHenson", "SimHensonRT", "SimGaussian"
type	NIGIC for the three Henson params for "SimHenson", "SimHensonRT". Type 'X' to specify your own A and B values (eg different dB scale)
A	parameter A for "SimHenson", "SimHensonRT"
B	parameter B for "SimHenson", "SimHensonRT"
cap	dB value for capping stdev form Henson formula for "SimHenson", "SimHensonRT"
rtData	data.frame with colnames == "Rt", "Dist", "Person" for "SimHensonRT"
rtFP	response time for false positives ??? for "SimHensonRT"

## Details

### Compass:

opiInitialize(ip, port)

If the chosen OPI implementation is Compass, then you must specify the IP address and port of the Compass server.

- ip is the IP address of the Compass server as a string.

- port is the TCP/IP port of the Compass server as a number. Warning: this returns a list, not a single error code.

**Display:** `opiInitialize((width, height, ppi, viewdist, lut = .OpiEnv$Display$LUT)`  
`)`

If the chosen OPI implementation is `Display`, then you can specify the limits of the plot area and the background color of the plot area. Note that this assumes `link{X11()}` is available on the platform.

We need to know the physical dimensions of the screen and the window generated in order to calculate stimulus position and size in degrees of visual angle. The physical dimensions in inches are calculated from width, height, and ppi. The pixel size pix per degree is then obtained using viewdist. A gamma function for the screen should be obtained and its lut passed to convert from luminance in cd/m<sup>2</sup> to 8-bit pixel value (256 levels).

## Value

Returns NULL if initialization succeeded, otherwise an implementation-dependent error.

### Compass:

Returns a list with elements: \* err NULL if successful, not otherwise. \* prl a pair giving the (x,y) in degrees of the Preferred Retinal Locus detected in the initial alignment. \* onh a pair giving the (x,y) in degrees of the ONH as selected by the user. \* image raw bytes being the JPEG compressed infra-red image acquired during alignment encoded using `openssl::base64_encode()`.

**Display:** Always returns NULL.

### Daydream:

Always returns NULL.

### imo:

Always returns NULL. Will stop if there is an error.

### Kowa AP-7000:

Always returns NULL.

### Octopus600:

Returns NULL if successful, or an Octopus 600 error code. The default background and stimulus setup is to white-on-white perimetry.

### Octopus900:

Returns NULL if successful, 1 if Octopus900 is already initialised by a previous call to `opiInitialize`, and 2 if some error occurred that prevented initialisation. The default background and stimulus setup is to white-on-white perimetry. Use `opiSetBackground` to change the background and stimulus colors.

### PhoneHMD:

Returns NULL if connection is made, otherwise, it returns a text with the error.

**Daydream**

```
opiInitialize(ip="127.0.0.1", port=50008, lut= seq(0, 400, length.out = 256), fovy = 90)
```

If the chosen OPI implementation is Daydream, then you must specify the IP address of the Android phone that is in the Daydream, and the port on which the server running on the phone is listening.

- ip is the IP address of the Daydream server as a string
- port is the TCP/IP port of the Daydream server as a number
- lut is a vector of 256 luminance values, with lut[i] being the cd/m<sup>2</sup> value for grey level i. Default is seq(0, 4000, length.out = 256)
- fovy Field of view in degrees in the y-axis. It is different depending on the device. For Daydream view, it is 90 degrees, for, Daydream view 2 is 100 degrees. Default is 90.

**imo**

```
opiInitialize(ip, port, ppp = 16, tracking = FALSE, tracktol = 2)
```

If the chosen OPI implementation is imo, then you must specify the IP address and port of the imo server.

- \* `{ip}` is the IP address of the imo server as a string.
- \* `{port}` is the TCP/IP port of the imo server as a number.
- \* `{ppp}` Pixel size in pixels per degree. Default is 16 ppp.
- \* `{tracking}` Whether to use tracking during stimulus presentation. Default is FALSE.
- \* `{tracktol}` Tolerance during tracking in degrees of visual angle. The system does not show any sti

**KowaAP7000**

```
opiInitialize(ip, port)
```

If the chosen OPI implementation is KowaAP7000, then you must specify the IP address and port of the AP-7000 server.

- ipAddress is the IP address of the AP-7000 server as a string.
- port is the TCP/IP port of the AP-7000 server as a number.

**Octopus600**

```
opiInitialize(ipAddress, eye, pupilTracking=FALSE, pulsar=FALSE, eyeControl=0)
```

If the chosen OPI implementation is Octopus600, then you must specify the IP address of the Octopus 600 and the eye to test.

ipAddress is the IP address of the Octopus 600 as a string.

eye must be either "left" or "right".

pupilTracking is TRUE to turn on IR illumination and set pupil black level (which happens at the first stimulus presentation).

pulsar is TRUE for pulsar stimulus, FALSE for size III white-on-white.

eyeControl \* 0 is off \* 1 is eye blink \* 2 is eye blink, forehead rest, fixation control \* 3 is eye blink, forehead rest, fixation control, fast eye movements

**Octopus900**

`opiInitialize(serverPort=50001,eyeSuiteSettingsLocation, eye, gazeFeed=NA, bigWheel=FALSE,pres_buzzer=0, resp_buzzer=0, zero_dB_is_10000_asb=TRUE)`

If the chosen OPI implementation is Octopus900, then you must specify a directory and the eye to be tested.

`serverPort` is the TCP/IP port on which the server is listening (on localhost).

`eyeSuiteSettingsLocation` is the folder name containing the EyeSuite setting files, and should include the trailing slash.

`eye` must be either "left" or "right".

`gazeFeed` is the name of an existing folder into which the video frames of eye tracker are recorded. Set to NA for no recording.

`bigWheel` is FALSE for a standard Octopus 900 machine. Some research machines are fitted with an alternate aperture wheel that has 24 sizes, which are accessed with `bigWheel` is TRUE. The mapping from size to 'hole on wheel' is hard coded; see code for details.

If `pres_buzzer` is greater than zero, a buzzer will sound with each stimuli presented.

If `resp_buzzer` is greater than zero, a buzzer will sound with each button press (response). The volume can be one of 0 (no buzzer), 1, 2, or 3 (max volume). If both buzzers are more than zero, the maximum of the two will be used as the volume.

If `zero_dB_is_10000_asb` is TRUE then 0 dB is taken as 10000 apostilbs, otherwise 0 dB is taken as 4000 apostilbs.

**PhoneHMD**

`opiInitialize(serverPort, port = 50008, lut = seq(0, 400, length.out = 256))` If the chosen OPI implementation is PhoneHMD, then you must specify the IP address of the Android PhoneHMD that is in the PhoneHMD, and the port on which the server running on the PhoneHMD is listening.

- `ip` is the IP address of the PhoneHMD server as a string
- `port` is the TCP/IP port of the PhoneHMD server as a number
- `lut` is a vector of 256 luminance values, with `lut[i]` being the  $\text{cd/m}^2$  value for grey level  $i$ . Default is `seq(0, 4000, length.out = 256)`

**SimGaussian**

`opiInitialize(sd, display=NA, maxStim=10000/pi)`

If the chosen OPI implementation is SimGaussian, then `sd` is the standard deviation value that the simulator will use for the slope/spread of the psychometric function.

`display` and `maxStim` is as for SimHenson.

**SimHenson**

`opiInitialize(type="C", A=NA, B=NA, cap=6, maxStim=10000/pi)`

If the chosen OPI implementation is SimHenson, then `type` can be one of: "N", for normal patients; "G", for POAG patients; and "C", for a combination. See Table 1 in Henson et al (2000).



If type is "X" then A and B should be specified and are used in place of one of the three A/B combinations as in Henson et al (2000). cap is the maximum standard deviation value that the simulator will use for the slope/spread of the psychometric function.

If display is a vector of four numbers c(xlow, xhi, ylow, yhi), then a plot area is created of dimension xlim=range(xlow, xhi) and ylim=range(ylow, yhi) and each call to opiPresent will display a point on the area. The color of the plot area can be set with opiSetBackground, and the color of the displayed point is determined by the stimulus passed to opiPresent.

maxStim is the maximum stimulus value in  $\text{cd/m}^2$ . This is used in converting  $\text{cd/m}^2$  to dB values, and vice versa.

### SimHensonRT

```
opiInitialize(type="C", A=NA, B=NA, cap=6, display=NA, maxStim=10000/pi, rtData, rtFP=1:1600)
```

If the chosen OPI implementation is SimHensonRT, then the first six parameters are as in SimHenson, and rtData is a data frame with at least 2 columns: "Rt", response time; and "Dist", signifying that distance between assumed threshold and stimulus value in your units.

This package contains RtSigmaUnits or RtDbUnits that can be loaded with the commands data(RtSigmaUnits) or data(RtDbUnits), and are suitable to pass as values for rtData.

rtFp gives the vector of values in milliseconds from which a response time for a false positive response is randomly sampled.

### References

David B. Henson, Shaila Chaudry, Paul H. Artes, E. Brian Faragher, and Alec Ansons. Response Variability in the Visual Field: Comparison of Optic Neuritis, Glaucoma, Ocular Hypertension, and Normal Eyes. Investigative Ophthalmology & Visual Science, February 2000, Vol. 41(2).

### See Also

[chooseOpi](#), [opiSetBackground](#), [opiClose](#), [opiPresent](#)

### Examples

```
## Not run:
# Set up the Compass
chooseOpi("Compass")
result <- opiInitialize(ip="192.168.1.7", port=44965)
if (is.null(result$err))
  print(result$prl)

## End(Not run)
## Not run:
# Set up a Display and wait for a key press in it.
chooseOpi("Display")
if (!is.null(opiInitialize(width = 1680, height = 1050, ppi = 128, viewdist = 25)))
  stop("opiInitialize failed")

opiSetBackground(lum = 100, color = "white", fixation = "Circle")
```

```

opiClose()

## End(Not run)
## Not run:
# Set up the imo
chooseOpi("imo")
opiInitialize(ip = "192.168.1.7", port = 1234)

## End(Not run)
## Not run:
# Set up the Kowa AP-7000
chooseOpi("KowaAP7000")
opiInitialize(ip="192.168.1.7", port=44965)

## End(Not run)
## Not run:
# Set up the Octopus 900
chooseOpi("Octopus900")
if (!is.null(opiInitialize(
  eyeSuiteSettingsLocation="C:/ProgramData/Haag-Streit/EyeSuite/",
  eye="left")))
  stop("opiInitialize failed")

## End(Not run)
# Set up a simulation using a psychometric function that is
# a cumulative gaussian of standard deviation 2
chooseOpi("SimGaussian")
if (!is.null(opiInitialize(sd=2)))
  stop("opiInitialize failed")
# Set up a simple simulation for white-on-white perimetry
chooseOpi("SimHenson")
if (!is.null(opiInitialize(type="C", cap=6)))
  stop("opiInitialize failed")

# Set up a simple simulation for white-on-white perimetry
# and display the stimuli in a plot region
chooseOpi("SimHenson")
if (!is.null(opiInitialize(type="C", cap=6)))
  stop("opiInitialize failed")
# Set up a simple simulation for white-on-white perimetry
# and display the stimuli in a plot region and simulate response times
chooseOpi("SimHensonRT")
data(RtSigmaUnits)
oi <- opiInitialize(type="C", cap=6, display=NA, rtData=RtSigmaUnits, rtFP=1:100)
if (!is.null(oi))
  stop("opiInitialize failed")

```

**Description**

List containing stimulus parameters with an S3 class attribute of opiKineticStimulus

**Usage**

"See details"

**Details**

The list should be of class opiKineticStimulus and contain the following elements.

- path list of (x,y) coordinates in degrees that is usable by xy.coords()
- image image[i] is the image to display (in a machine specific format) in the section of the path specified by path[i]..path[i+1].
- levels if is.na(image) then levels[i] is the stimulus level in cd/m<sup>2</sup> in the section of the path specified by path[i]..path[i+1]
- sizes sizes[i] is the size of stimulus (diameter in degrees) to use for the section of path specified by path[i]..path[i+1], or a scaling factor for images[i].
- colors colors[i] is the color to use for the stimulus in the section of path specified by path[i]..path[i+1]. Ignored if !is.na(image).
- speeds speeds[i] is the speed (degrees per second) for the stimulus to traverse the path specified by path[i]..path[i+1].
- ... machine specific parameters

**Octopus 900**

x and y are in degrees, with precision to three decimal places recognised.

image is not possible on an Octopus 900.

levels are in cd/m<sup>2</sup>, and are rounded to the nearest one tenth of a dB for display.

colors are ignored. Use opiSetBackground() to alter stimulus color.

sizes are in degrees, but are rounded to the nearest Goldmann Size I.V for display.

**Kowa AP 7000**

Only a simple path with a start and an end point is supported by the AP-7000.

x and y are in degrees and should only be length 2. (precision?)

image is not possible on an Kowa AP 7000.

levels are in cd/m<sup>2</sup> in the range 0.03 to 3183, and are rounded to the nearest one tenth of a dB for display. (precision?)

colors one of .OpiEnv\$KowaAP7000\$COLOR\_WHITE, .OpiEnv\$KowaAP7000\$COLOR\_GREEN, .OpiEnv\$KowaAP7000\$COLOR\_ and .OpiEnv\$KowaAP7000\$COLOR\_RED.

sizes are in degrees, but are rounded to the nearest Goldmann Size I.V for display.

speeds are in degrees per second in the range 3 to 5.

**Compass**

Not implemented.

**See Also**

[opiSetBackground](#), [opiStaticStimulus](#), [opiTemporalStimulus](#)

**Examples**

```
# A Size III white kinetic stimuli on a bilinear path {(27,27), (15,20), (0,0)}
stim <- list(path=list(x=c(27,15,0), y=c(27,20,0)),
            izes=rep(0.43,2),
            colors=rep("white",2),
            levels=rep(318,2),
            speeds=c(4,3))
class(stim) <- "opiKineticStimulus"
```

---

 opiPresent

*Use OPI to present stimulus*


---

**Description**

Generic function for presentation of stimulus stim. Depending on your choice of OPI implementation set using `chooseOpi()`, different parameters are available for `opiPresent`

**Usage**

```
opiPresent(stim, nextStim = NULL, ...)

compass.opiPresent(stim, nextStim = NULL)

display.opiPresent(stim, nextStim = NULL)

daydream.opiPresent(stim, nextStim = NULL)

imo.opiPresent(stim, nextStim = NULL)

kowaAP7000.opiPresent(stim, nextStim = NULL)

octo600.opiPresent(stim, nextStim = NULL)

octo900.opiPresentF310(stim, nextStim = NULL)

PhoneHMD.opiPresent(stim, nextStim = NULL)

simG.opiPresent(stim, nextStim = NULL, fpr = 0.03, fnr = 0.01, tt = 30)
```

```

simH.opiPresent(
    stim,
    nextStim = NULL,
    fpr = 0.03,
    fnr = 0.01,
    tt = 30,
    criteria = 0.97,
    rt_shape = 5.3,
    rt_rate = 1.4,
    rt_scale = 0.1
)

simH_RT.opiPresent(
    stim,
    nextStim = NULL,
    fpr = 0.03,
    fnr = 0.01,
    tt = 30,
    notSeenToSeen = TRUE
)

simNo.opiPresent(stim, nextStim = NULL)

simYes.opiPresent(stim, nextStim = NULL)

```

### Arguments

stim	a list of class <a href="#">opiStaticStimulus</a> , <a href="#">opiKineticStimulus</a> , or <a href="#">opiTemporalStimulus</a> to be presented.
nextStim	unused - included for compliance with OPI standard.
...	Parameters specific to your chosen opi implementation
fpr	false positive rate for OPI implementation "SimHenson"
fnr	false negative rate for OPI implementation "SimHenson"
tt	SOMETHING for OPI implementation "SimHenson"
criteria	CRITERIA for OPI implementation "SimHenson"
rt_shape	response time shape parameter for OPI implementation "SimHenson"
rt_rate	response time rate parameter for OPI implementation "SimHenson"
rt_scale	response time scale parameter for OPI implementation "SimHenson"
notSeenToSeen	SOMETHING for OPI implementation "SimHensonRT"

### Details

[opiPresent](#) is blocking in that it will not return until either a response is obtained, or at least the responseWindow milliseconds has expired. (Note that more time might have expired.) Specifying nextStim allows the implementing machine to use the time waiting for a response to stim to make preparations for the next stimuli. (For example retargeting the projector or moving aperture and/or

filter wheels.) There is no guarantee that the next call to `opiPresent` will have `nextStim` as the first argument; this could be checked by the machine specific implementations (but currently is not, I think).

Also note that to allow for different parameters depending on the implementation chosen with `chooseOpi`, every parameter MUST be named in a call to `opiPresent`.

**Display:** Present a circle of radius `stim$size` and color `stim$color` at (`stim$x`, `stim$y`) for `stim$duration` ms and wait for a keyboard or mouse response for `stim$responseWindow` ms. `stim$size`, `stim$x` and `stim$y` are in the same units as `xlim` and `ylim` as specified in `opiInitialise`. If the chosen OPI implementation is `Display`, then `nextStim` is ignored. Duration and response window are rounded to the nearest 5 ms. Currently only implemented for `opiStaticStimulus`.

**Daydream:** If the chosen OPI implementation is `Daydream`, then `nextStim` is ignored. Note that the dB level is rounded to the nearest  $\text{cd/m}^2$  that is in the lut specified in `opiInitialise`. Currently uses the most simple algorithm for drawing a 'circle' (ie not Bresenham's). Currently only implemented for `opiStaticStimulus`.

**imo:** DETAILS HERE

**KowaAP7000:** `opiPresent(stim, nextStim=NULL)`

If the chosen OPI implementation is `KowaAP7000`, then `nextStim` is ignored.

## Value

A list containing

- `err` NULL if no error occurred, otherwise a machine-specific error message. This should include errors when the specified size cannot be achieved by the device (for example, in a projection system with an aperture wheel of predefined sizes.) If `stim` is NULL, then `err` contains the status of the machine.
- `seen` TRUE if a response was detected in the allowed `responseWindow`, FALSE otherwise. (Note, see `Octopus900F310` above).
- `time` The time in milliseconds from the onset (or offset, machine-specific) of the presentation until the response from the subject if `seen` is TRUE. If `seen` is FALSE, this value is undefined. For kinetic perimetry on the `O900`, this value is unknown... (what does this mean!?)

## O600:

answer only returned for `Octopus600`. Can be the following values:

- 0 = stimulus not seen;
- 1 = stimulus seen;
- 132 = Response button was pressed before stimulus presentation (Patient needs a break - hold on examination);
- 36 = Eye is closed before stimulus presentation;
- 68 = Fixation lost before stimulus presentation (pupil center is out of green window in video image);
- 260 = Forehead rest lost before stimulus presentation;

- 516 = Fast Eye movements before stimulus presentation;
- 258 = Forehead rest lost during stimulus presentation;
- 66 = Fixation lost during stimulus presentation (pupil center is out of green window in video image);
- 34 = Eye was closed during stimulus presentation;
- 18 = Patient answer was too early ( $\leq 100$ ms after stimulus presentation) - lucky punch;
- 514 = Fast Eye movements during stimulus presentation

#### **Kowa AP7000:**

- `pupilX` Only returned for KowaAP7000 (in pixels) and an `opiStaticStimulus` or O900 (in degrees) and static/kinetic if `gazeFeed==1`. x-coordinate of centre of pupil during presentation.
- `pupilY` Only returned for KowaAP7000 (in pixels) and an `opiStaticStimulus` or O900 (in degrees) and static/kinetic if `gazeFeed==1`. y-coordinate of centre of pupil during presentation.
- `purkinjeX` Only returned for KowaAP7000 and an `opiStaticStimulus`. x-coordinate of centre of Purkinje Image in pixels during presentation.
- `purkinjeY` Only returned for KowaAP7000 and an `opiStaticStimulus`. y-coordinate of centre of Purkinje Image in pixels during presentation.

#### **Kowa AP7000 and Octopus O900:**

- `x` Only returned for KowaAP7000 or Octopus900 and an `opiKineticStimulus`. x coordinate of stimuli when button is pressed.
- `y` Only returned for KowaAP7000 or Octopus900 and an `opiKineticStimulus`. y coordinate of stimuli when button is pressed.

#### **Compass:**

- `time_rec` Time since epoch that the `opiPresent` command was received by the Compass in ms.
- `time_hw` Hardware time of button press or response window expired (integer ms). To get the hardware time that a presentation began, subtract `responseWindow` from `th` (for aligning with fixation data returned by `opiClose()`).
- `time_resp` Time since epoch that the response was received or response window expired (in ms).
- `num_track_events` The number of tracking events associated with this presentation.
- `num_motor_fails` The number of time the motor could not keep pace with eye movements.
- `pupil_diam` The diameter of the pupil on millimeters on presentation.
- `loc_x` The x location in pixels of the presentation on the retinal image returned by `opiInitialize`.
- `loc_y` The y location in pixels of the presentation on the retinal image returned by `opiInitialize`.

#### **Compass:**

A list containing

- `err` 0 all clear,  $\geq 1$  some error codes (eg cannot track, etc) (integer)
- `seen` FALSE for not seen, TRUE for seen (button pressed in response window)
- `time` response time in ms (integer) since stimulus onset, -1 for not seen
- `time_rec` time since epoch when command was received at Compass (integer ms)
- `time_pres` time since epoch that stimulus was presented (integer ms)

- num\_track\_events number of tracking events that occurred during presentation (integer)
- num\_motor\_fails number of times motor could not follow fixation movement during presentation (integer)
- pupil\_diam pupil diameter in mm (float)
- loc\_x pixels integer, location in image of presentation (integer)
- loc\_y pixels integer, location in image of presentation (integer)

### Compass

opiPresent(stim, nextStim=NULL)

If the chosen OPI implementation is Compass, then nextStim is ignored. Note that the dB level is rounded to the nearest integer.

If tracking is on, then this will block until the tracking is obtained, and the stimulus presented.

### Octopus600

opiPresent(stim, nextStim=NULL)

If the chosen OPI implementation is Octopus600, then nextStim is ignored. If eyeControl is non-zero, as set in opiInitialize, answer codes describing patient state may arise (see answer field in the Value section).

### Octopus900F310

opiPresent(stim, nextStim=NULL)

This functions as for the Octopus900, but responses are taken from the F310 Controller.

If the L button is pressed, seen is set to 1.

If the R button is pressed, seen is set to 2.

If no button is pressed within responseWindow, then seen is set to 0.

### PhoneHMD

If the chosen OPI implementation is PhoneHMD, then nextStim is ignored. PhonVR

### SimGaussian

opiPresent(stim, nextStim=NULL, fpr=0.03, fnr=0.01, tt=30)

If the chosen OPI implementation is SimGaussian, then the response to a stimuli is determined by sampling from a Frequency-of-Seeing (FoS) curve (also known as the psychometric function) with formula  $fpr + (1 - fpr - fnr) * (1 - pnorm(x, tt, simG.global.sd))$ , where  $x$  is the stimulus value in Humphrey dB, and  $simG.global.sd$  is set with opiInitialize.



**SimHenson**

```
opiPresent(stim, nextStim=NULL, fpr=0.03, fnr=0.01, tt=30)
```

If the chosen OPI implementation is SimHenson, then the response to a stimuli is determined by sampling from a Frequency-of-Seeing (FoS) curve (also known as the psychometric function) with formula

$$fpr + (1 - fpr - fnr)(1 - pnorm(x, tt$$

, where  $x$  is the stimulus value in Humphrey dB, and  $pxVar$  is

$$\min\left(\text{simH.global.cap}, e^{A \times tt + B}\right).$$

The ceiling `simH.global.cap` is set with the call to `opiInitialize`, and  $A$  and  $B$  are from Table 1 in Henson et al (2000). Which values are used is determined by `simH.type` which is also set in the call to `opiInitialize`.

Note that if the stimulus value is less than zero, then the Henson formula is not used. The probability of seeing is `fpr`.

```
opiPresent(stim, nextStim=NULL, fpr=0.03, fnr=0.01, tt=NULL, criteria=0.95, rt_shape=5.3,
rt_rate=1.4, rt_scale=0.1)
```

For determining seen/not-seen for kinetic, the first location (to a fidelity of 0.01 degrees) on the path (it only works for single paths now) where the probability of seeing is equal to `criteria` is found. If no such location exists, then the stimuli is not seen. The probability of seeing at each location is determined using a frequency-of-seeing curve defined as a cumulative Gaussian with parameters controlled by `tt` and `opiInitialize`. At each location along the path, the mean of the FoS is taken from the `tt` function, which takes a distance-along-path (in degrees) as an argument, and returns a dB value which is the static threshold at that distance along the path.

Function `tt` can return NA for not thresholds that are always not seen. At each location along the path, the standard deviation of the FoS is sampled from a Gaussian with mean taken from the formula of Henson et al (2000), as parametrised by `opiInitialize`, and standard deviation 0.25.

The location of a false positive response (for the total kinetic path) is sampled uniformly from the start of the path to the 'seeing' location, or the entire path if the stimuli is not seen.

Note that the false positive rate `fpr` and the false negative rate `fnr` are specified for the whole path, and not for the individual static responses along the way.

The actual location returned for a seen response is the location where the probability of seeing equals `criteria`, plus a response time sampled from a Gamma distribution parameterised by `rt_shape` and `rt_rate` and multiplied by `rt_scale`. That is: `rgamma(1, shape=rt_shape, rate=rt_rate) / rt_scale`.

**SimHensonRT**

```
opiPresent(stim, nextStim=NULL, fpr=0.03, fnr=0.01, tt=30, dist=stim$level - tt)
```

For static stimuli, this function is the same as for SimHenson, but reaction times are determined by sampling from `rtData` as passed to `opiInitialize`. The `dist` parameter is the distance of the stimulus level from the true threshold, and should be in the same units as the `Dist` column of `rtData`. The default is just the straight difference between the stimulus level and the true threshold, but you might want it scaled somehow to match `rtData`.

**SimNo**

```
opiPresent(stim, nextStim=NULL)
```

If the chosen OPI implementation is SimNo, then the response to a stimuli is always no, hence [opiPresent](#) always returns `err=NULL`, `seen=FALSE`, and `time=0`.

**SimYes**

```
opiPresent(stim, nextStim=NULL)
```

If the chosen OPI implementation is SimYes, then the response to a stimuli is always yes, hence [opiPresent](#) always returns `err=NULL`, `seen=TRUE`, and `time=0`.

**References**

David B. Henson, Shaila Chaudry, Paul H. Artes, E. Brian Faragher, and Alec Ansons. Response Variability in the Visual Field: Comparison of Optic Neuritis, Glaucoma, Ocular Hypertension, and Normal Eyes. *Investigative Ophthalmology & Visual Science*, February 2000, Vol. 41(2).

**See Also**

[opiStaticStimulus](#), [opiKineticStimulus](#), [opiTemporalStimulus](#), [chooseOpi](#), [opiInitialize](#)

**Examples**

```
## Not run:
# Display a spot
chooseOpi("Display")
if(!is.null(opiInitialize(width = 1680, height = 1050, ppi = 128, viewdist = 25)))
  stop("opiInitialize failed")
opiSetBackground(lum = 50, color = "white", fixation = "Cross")

makeStim <- function(db) {
  s <- list(x = 9, y = 9, level = dbTocd(db, 400), size = 1.72, color = "white",
           duration = 1000, responseWindow = 1000)
  class(s) <- "opiStaticStimulus"
  return(s)
}
result <- opiPresent(makeStim(0))

opiClose()

## End(Not run)
# Stimulus is Size III white-on-white as in the HFA
makeStim <- function(db, n) {
  s <- list(x=9, y=9, level=dbTocd(db, 10000/pi), size=0.43, color="white",
           duration=200, responseWindow=1500)
  class(s) <- "opiStaticStimulus"
  return(s)
}

chooseOpi("SimHenson")
```

```
if (!is.null(opiInitialize(type="C", cap=6)))
  stop("opiInitialize failed")

result <- opiPresent(stim=makeStim(10,0), tt=30, fpr=0.15, fnr=0.01)

# Will not work as 'stim' is not named
#result <- opiPresent(makeStim(10,0), tt=30, fpr=0.15, fnr=0.01)

if (!is.null(opiClose()))
  warning("opiClose() failed")
# Same but with simulated reaction times
chooseOpi("SimHensonRT")
data(RtSigmaUnits)
if (!is.null(opiInitialize(type="C", cap=6, rtData=RtSigmaUnits)))
  stop("opiInitialize failed")

dist <- (10 - 30)/min(exp(-0.098 * 30 + 3.62), 6)
result <- opiPresent(stim=makeStim(10,0), tt=30, fpr=0.15, fnr=0.01, dist=dist)

if (!is.null(opiClose()))
  warning("opiClose() failed")
```

---

opiQueryDevice

*Query device using OPI*

---

## Description

Generic function for getting details of the chosen OPI implementation that is set with `chooseOpi()`

## Usage

```
opiQueryDevice(...)

compass.opiQueryDevice()

display.opiQueryDevice()

daydream.opiQueryDevice()

imo.opiQueryDevice()

kowaAP7000.opiQueryDevice()

octo600.opiQueryDevice()

octo900.opiQueryDevice()

PhoneHMD.opiQueryDevice()

simG.opiQueryDevice()
```

**Arguments**

... Implementation specific parameters. See details.

**Details**

**Compass:** Return a list of all the constants used in the OPI Compass module.

**Display:** Returns all constants in `.OpiEnv$Display` as a list.

**Display:** Returns values in use by Display.

**Daydream:** Returns all constants in `.OpiEnv$DayDream` as a list.

**Daydream:** DETAILS

**Value**

Returns a list that contains `isSim` and implementation-dependent data.

`isSim` is TRUE if the device is a simulation, or FALSE if the device is a physical machine.

**Compass:** A list containing constants and their value used in the OPI Compass module.

**Octopus600:**

Returns a list of 10 items: \* `answerButton` {0 = not pressed, 1 = pressed } \* `headSensor` {0 = no forehead detected, 1 = forehead detected } \* `eyeLidClosureLeft` {0 = eye is open, 1 = eye is closed} \* `eyeLidClosureRight` {0 = eye is open, 1 = eye is closed} \* `fixationLostLeft` {1 = eye pos lost, 0 = eye pos ok} \* `fixationLostRight` {1 = eye pos lost, 0 = eye pos ok} \* `pupilPositionXLeft` (in px) \* `pupilPositionYLeft` (in px) \* `pupilPositionXRight` (in px) \* `pupilPositionYRight` (in px)

**Octopus900:**

list containing `isSim=FALSE`.

**KowaAP7000**

If the chosen OPI is `KowaAP7000`, then this function returns the current location of the pupil. See the Value section for details.

Returns a list of 4 items:

- `pupilX`, the x-coordinate of the pupil position in pixels.
- `pupilY`, the y-coordinate of the pupil position in pixels.
- `purkinjeX`, the x-coordinate of the purkinje position in pixels.
- `purkinjeY`, the y-coordinate of the purkinje position in pixels.

It also prints a list of constants that OPI knows about for the AP-7000.

**Octopus600**

If the chosen OPI is `Octopus600`, then this function returns information about the patient. See the Value section for details.

**Octopus900**

Prints defined constants in OPI package pertaining to Octopus 900.

**PhoneHMD**

Returns all constants in `.OpiEnv$PhoneHMD` as a list.

**See Also**

[chooseOpi](#)

**Examples**

```
chooseOpi("SimGaussian")
if (!is.null(opiInitialize(sd=2)))
  stop("opiInitialize failed")
print(opiQueryDevice())
```

---

<code>opiSetBackground</code>	<i>Set background using OPI</i>
-------------------------------	---------------------------------

---

**Description**

Generic function for setting background of the chosen OPI implementation that is set with `chooseOpi()`

**Usage**

```
opiSetBackground(...)

compass.opiSetBackground(lum = NA, color = NA, fixation = NA, tracking_on = NA)

display.opiSetBackground(
  lum = .OpiEnv$Display$background_lum,
  color = .OpiEnv$Display$background_color,
  fixation = .OpiEnv$Display$fixation,
  fix_cx = .OpiEnv$Display$fix_cx,
  fix_cy = .OpiEnv$Display$fix_cy,
  fix_sx = .OpiEnv$Display$fix_sx,
  fix_sy = .OpiEnv$Display$fix_sy,
  fix_color = .OpiEnv$Display$fix_color
)

daydream.opiSetBackground(
  eye,
  lum = 10,
  color = "white",
  fixation = "None",
```

```
    fix_cx = 0,  
    fix_cy = 0,  
    fix_sx = 2,  
    fix_sy = 2,  
    fix_color = "green"  
  )  
  
imo.opiSetBackground(  
  bgl = 30,  
  tgl = 1,  
  tgl db = 25,  
  tglx = 0,  
  tgly = 0,  
  bgr = 30,  
  tgr = 1,  
  tgr db = 25,  
  tgrx = 0,  
  tgry = 0  
)  
  
kowaAP7000.opiSetBackground(lum = NA, color = NA, fixation = NA)  
  
octo600.opiSetBackground(  
  bgColor = NA,  
  fixType = NA,  
  fixColor = NA,  
  fixIntensity = 255  
)  
  
octo900.opiSetBackground(  
  lum = NA,  
  color = NA,  
  fixation = NA,  
  fixIntensity = NA  
)  
  
PhoneHMD.opiSetBackground(  
  bgeye,  
  bglum = 10,  
  bgcol = "white",  
  fixeye,  
  fixtype = "none",  
  fixcx = 0,  
  fixcy = 0,  
  fixsx = 2,  
  fixsy = 2,  
  fixtheta = 0,  
  fixlum = 100,
```

```

    fixcol = "green"
)

simG.opiSetBackground(col, gridCol)

simH.opiSetBackground(col, gridCol)

simH_RT.opiSetBackground(col, gridCol)

simNo.opiSetBackground(col, gridCol)

simYes.opiSetBackground(col, gridCol)

```

### Arguments

...	Implementation specific parameters. See details.
lum	Luminance level in cd/m <sup>2</sup>
color	Stimulus color
fixation	fixation target
tracking_on	TRUE for tracking on, FALSE for off
fix_cx	fixation x position in degrees of visual angle. Default is 0
fix_cy	fixation y position in degrees of visual angle. Default is 0
fix_sx	fixation horizontal size in degrees of visual angle. Default is 1
fix_sy	fixation vertical size in degrees of visual angle. Default is 1
fix_color	fixation color
eye	eye
bgl	left eye background luminance in dB Default is 25
tgl	left eye fixation target. Default is 1
tgldb	left eye fixation target luminance in dB. Default is 20
tg1x	left eye fixation target x-position in degrees. Default is 0.
tg1y	left eye fixation target y-position in degrees. Default is 0.
bgr	right eye background luminance in dB. Default is 25
tgr	right eye target type. Default is 1
tgrdb	right eye target luminance in dB. Default is 20
tgrx	right eye fixation target x-position in degrees. Default is 0.
tgry	right eye fixation target y-position in degrees. Default is 0.
bgColor	Background color
fixType	fixation type
fixColor	fixation color
fixIntensity	fixation point intensity
bgeye	eye where to display the background, can be left, right, or both

bglum	background luminance
bgcol	background color. Either a color name that R can understand or a vector with R, G, B, and alpha channels
fixeye	eye for the fixation target, can be left, right, or both
fixtype	the target of the fixation
fixcx	fixation x position in degrees of visual angle. Default is 0
fixcy	fixation y position in degrees of visual angle. Default is 0
fixsx	fixation horizontal size in degrees of visual angle. Default is 1
fixsy	fixation vertical size in degrees of visual angle. Default is 1
fixtheta	angle of rotation of the fixation target in degrees. Default is 0
fixlum	fixation luminance
fixcol	fixation color. Same specifications as background color
col	DESCRIPTION for OPI implementations based on simulations
gridCol	DESCRIPTION for OPI implementation based on simulations

### Details

**Compass:** `opiSetBackground(fixation=NA, tracking_on=NA)`

- `fixation=c(x,y,t)` where
  - x is one of -20, -6, -3, 0, 3, 6, 20 degrees.
  - y is 0 degrees.
  - t is 0 for a spot fixation marker at  $c(x, y)$ , or 1 for a square centred on one of  $(-3, 0)$ ,  $(0, 0)$ ,  $(+3, 0)$ .
- `tracking_on` is either 0 (tracking off) or 1 (tracking on).

Note: tracking will be relative to the PRL established with the fixation marker used at setup (call to OPI-OPEN), so when tracking is on you should use the same fixation location as in the setup.

**imo:** DETAILS

### Value

Returns NULL if succeeded, otherwise an implementation-dependent error as follows.

**Compass:** A list containing error which is NULL for success, or some string description for fail.

**Display:**

Changes the background and the fixation marker.

**Daydream:**

Returns NULL or an error string.

**imo:** DETAILS

**KowaAP7000:**

Always returns NULL



**Octopus600:**

- -1 to be implemented
- -2 O600 sent back an error; bad background parameters
- -3 O600 sent back an error; bad fixation parameters
- NULL Success

**Octopus900:**

- -1 indicates opiInitialize has not been called.
- -2 indicates could not set the background color.
- -3 indicates could not set the fixation marker.
- -4 indicates that all input parameters were NA.

**PhoneHMD:**

Sets background for left, right, or both eyes in PhoneHMD.

**Display**

opiSetBackground(TODO)

**Daydream**

opiSetBackground(eye, lum=10, color="white", fixation="Cross", fix\_cx=0, fix\_cy=0, fix\_sx=2, fix\_sy=2, fix\_lum=10, fix\_color="green")

- lum background luminance in  $\text{cd/m}^2$  is set to nearest grey value in lut from opiInitialize. Default is  $10 \text{ cd/m}^2$
- color color of the background. It can be 'white' (default) or 'green'.
- fixation can only be 'Cross' at the moment.
- fix\_cx, fix\_cy fixation (x, y) position in degrees of visual angle.
- fix\_sx, fix\_sy dimensions of fixation target in degrees of visual angle.
- fix\_lum luminance of the fixation target in  $\text{cd/m}^2$  is set to nearest grey value in lut from opiInitialize. Default is  $15 \text{ cd/m}^2$ .
- fix\_color color of the fixation target. It can be 'white' or 'green' (default).

**KowaAP7000**

opiSetBackground(lum, color, fixation)

lum and color are dependant for the Kowa AP-7000. A white background must be  $10 \text{ cd/m}^2$ , and a yellow background must be  $100 \text{ cd/m}^2$ .

If lum is 10 and color is not set, then .OpiEnv\$KowaAP7000\$BACKGROUND\_WHITE is assumed.

If lum is 100 and color is not set, then .OpiEnv\$KowaAP7000\$BACKGROUND\_YELLOW is assumed.

If both lum and color is set, then lum is ignored (a warning will be generated if lum is incompatible with color).

fixation is one of \* .OpiEnv\$KowaAP7000\$FIX\_CENTER, fixation marker in the centre. \* .OpiEnv\$KowaAP7000\$FIX\_CENTR fixation marker in the centre. \* .OpiEnv\$KowaAP7000\$FIX\_AUX, fixation marker is ????. \* .OpiEnv\$KowaAP7000\$FIX\_MACUL fixation marker is a circle(?). \* .OpiEnv\$KowaAP7000\$FIX\_AUX\_LEFT, fixation marker is as for AUX but only lower left.

**Octopus600**

This function has no effect.

**Octopus900**

`opiSetBackground(lum=NA, color=NA, fixation=NA, fixIntensity=NA)`

`lum` is intensity of the background and can be one of \* `.OpiEnv$0900$BG_OFF`, which turns background off. \* `.OpiEnv$0900$BG_1`, background of 1.27 cd/m<sup>2</sup>. \* `.OpiEnv$0900$BG_10`, background of 10 cd/m<sup>2</sup>. \* `.OpiEnv$0900$BG_100`, background of 100 cd/m<sup>2</sup>.

`color` can be one of the following choices. \* `.OpiEnv$0900$MET_COL_WW` for white-on-white \* `.OpiEnv$0900$MET_COL_RW` for red-on-white \* `.OpiEnv$0900$MET_COL_BW` for blue-on-white \* `.OpiEnv$0900$MET_COL_WY` for white-on-yellow \* `.OpiEnv$0900$MET_COL_RY` for red-on-yellow \* `.OpiEnv$0900$MET_COL_BY` for blue-on-yellow

`fixation` is one of \* `.OpiEnv$0900$FIX_CENTRE` or `.OpiEnv$0900$FIX_CENTER` \* `.OpiEnv$0900$FIX_CROSS` \* `.OpiEnv$0900$FIX_RING`

`\code{fixIntensity}` is a percentage between 0 and 100. 0 is off, 100 the brightest.

Note if you specify fixation you also have to specify `fixIntensity`.

**PhoneHMD**

`opiSetBackground(bgeye, bglum = 10, bgcol = "white", fixeye = eye, fixtype = "Cross", fixlum = 100, fixcol = "green", fixcx = 0, fixcy = 0, fixsx = 2, fixsy = 2)`

- `bgeye` eye for the background. Can be "R", "L", or "B" for right, left, or both
- `bglum` background luminance in cd/m<sup>2</sup> is set to nearest grey value in lut from `opiInitialize`. Default is 10 cd/m<sup>2</sup>
- `bgcol` color of the background. Default is white
- `fixeye` eye for the background. Takes the same values as `eye` (default)
- `fixtype` fixation target
- `fixcx`, `fixcy` fixation (x, y) position in degrees of visual angle
- `fixsx`, `fixsy` dimensions of fixation target in degrees of visual angle
- `fixtheta` angle of rotation of the fixation target
- `fixlum` luminance of the fixation target. Default is 100 cd/m<sup>2</sup>
- `fixcol` color of the fixation target. Default is green

**SimGaussian**

`opiSetBackground(col, gridCol)`

`col` is the background color of the plot area used for displaying stimuli, and `gridCol` the color of the gridlines. Note the plot area will only be displayed if `opiInitialize` is called with a valid `display` argument.

**SimHenson**

```
opiSetBackground(col, gridCol)
```

col is the background color of the plot area used for displaying stimuli, and gridCol the color of the gridlines. Note the plot area will only be displayed if opiInitialize is called with a valid display argument.

**SimHensonRT**

```
opiSetBackground(col, gridCol)
```

col is the background color of the plot area used for displaying stimuli, and gridCol the color of the gridlines. Note the plot area will only be displayed if opiInitialize is called with a valid display argument.

**SimNo**

DETAILS

**SimYes**

DETAILS

**See Also**

[chooseOpi](#)

**Examples**

```
## Not run:
# Set up a Display and wait for a key press in it.
chooseOpi("Display")
if (!is.null(opiInitialize(width = 1680, height = 1050, ppi = 128, viewdist = 25)))
  stop("opiInitialize failed")

opiSetBackground()
opiSetBackground(lum = 100, color = "white", fixation = "Circle")
opiSetBackground(lum = 100, color = "white", fixation = "Cross", fix_color = "red",
                  fix_cx = 2, fix_cy = 5, fix_sx = 1, fix_sy = 2)
opiClose()

## End(Not run)
## Not run:
chooseOpi("Octopus900")
oi <- opiInitialize(eyeSuiteJarLocation="c:/EyeSuite/",
                   eyeSuiteSettingsLocation="c:/Documents and Settings/All Users/Haag-Streit/",
                   eye="left")
if(!is.null(oi))
  stop("opiInitialize failed")
if(!is.null(opiSetBackground(fixation=.OpiEnv$0900$FIX_CENTRE)))
  stop("opiSetBackground failed")
if(!is.null(opiSetBackground(fixation=.OpiEnv$0900$FIX_RING, fixIntensity=0)))
```

```

    stop("opiSetBackground failed")
  if(!is.null(opiSetBackground(color=.OpiEnv$0900$MET_COL_BY)))
    stop("opiSetBackground failed")
  if(!is.null(opiSetBackground(lum=.OpiEnv$0900$BG_100, color=.OpiEnv$0900$MET_COL_RW)))
    stop("opiSetBackground failed")
  opiClose()

## End(Not run)
chooseOpi("SimGaussian")
if (!is.null(opiInitialize(sd=2)))
  stop("opiInitialize failed")
if (!is.null(opiSetBackground(col="white",gridCol="grey")))
  stop("opiSetBackground failed, which is very surprising!")

```

---

opiStaticStimulus      *Stimulus parameter list*

---

### Description

List containing stimulus parameters with an S3 class attribute of opiStaticStimulus

### Usage

"See details"

### Details

The list should be of class opiStaticStimulus and contain the following elements.

- x coordinate of the center of stimulus in degrees relative to fixation
- y coordinate of the center of stimulus in degrees relative to fixation
- image an image to display in a machine specific format
- level stimulus level in  $\text{cd/m}^2$  (ignored if !is.na(image))
- size diameter of target in degrees, or scaling factor for image if specified
- color machine specific stimulus color settings (ignored if !is.na(image))
- duration total stimulus duration in milliseconds maximum responseWindow time ( $\geq 0$ ) in milliseconds to wait for a response from the onset of the stimulus presentation
- ... machine-specific parameters

### SimHenson and SimGaussian

Only level is used.

Duration and location are ignored, color is assumed "white" and size is assumed to be 26/60 (Goldmann III).

**Octopus 900**

x and y are in degrees, with precision to one decimal place recognised.

image is not possible on an Octopus 900.

level is in  $\text{cd/m}^2$ , and is rounded to the nearest one tenth of a dB for display.

color is ignored. Use `opiSetBackground()` to alter stimulus color

`checkFixationOK` is a function that takes the return value from `opiPresent` and returns either TRUE, indicating that fixation was good for the presentation; or FALSE, indicating that fixation was not good for the presentation.

**Octopus 900 F310 Controller**

As for the Octopus 900, but a `responseWindow` of -1 means that the Octopus 900 server will wait until either the L and R button is pressed in the controller until returning.

**Kowa AP7000**

x and y are in degrees. (precision?)

image is not possible on an Kowa AP 7000.

level are in  $\text{cd/m}^2$  in the range 0.03 to 3183, nearest one tenth of a dB for display.

size is in degrees, but is rounded to the nearest Goldmann Size I.V for display.

color one of `.OpiEnv$KowaAP7000$COLOR_WHITE`, `.OpiEnv$KowaAP7000$COLOR_GREEN`, `.OpiEnv$KowaAP7000$COLOR_BLUE` and `.OpiEnv$KowaAP7000$COLOR_RED`.

**imo**

x, y, level, size, and color are not used.

image is a list of two matrices: the first for the right eye, the second for the left. Each image is a 1080x1080 matrix with each element in the range 0 to 80, which maps onto 0dB to 40dB in steps of 0.5dB. Thus 0 is 0dB, 3283.048  $\text{cd/m}^2$ ; 1 is 0.5dB; and 80 is 40dB, 10  $\text{cd/m}^2$ .

tracking is TRUE if auto image placement to keep pupil centred is used, or FALSE to turn off imo auto-image placement to keep centred on pupil.

**Compass**

x and y are in degrees (floating point) (range -30 to 30 inclusive).

level is in  $\text{cd/m}^2$ , and is rounded to the nearest whole dB for display (range 0 to 50). 0dB is 10000aps.

`responseWindow` is in milliseconds (range 0 to 2680). Parameter duration is assumed to be 200ms, size is assumed to be Goldmann III (0.43), and color is assumed to be white.

**See Also**

[opiSetBackground](#), [opiKineticStimulus](#), [opiTemporalStimulus](#)

**Examples**

```
stim <- list(x=9, y=9, image=NA, 314, size=0.43, color="white",
            duration=200, responseWindow=1500)
class(stim) <- "opiStaticStimulus"
```

---

opiTemporalStimulus     *Stimulus parameter list*

---

**Description**

List containing stimulus parameters with an S3 class attribute of opiTemporalStimulus

**Usage**

"See details"

**Details**

The list should be of class opiTemporalStimulus and contain the following elements.

- x coordinate of the center of stimulus in degrees relative to fixation
- y coordinate of the center of stimulus in degrees relative to fixation
- image an image to display in a machine specific format
- lut if `is.na(image)` then this is a lookup table (vector) for stimulus level at each step of rate Hz in  $\text{cd/m}^2$ . If image is specified, then this is a list of images, in the same format as image, that is stepped through at rate Hz.
- size diameter of target in degrees, or scaling factor for image if specified
- color machine specific stimulus color settings (ignored if `!is.na(image)`)
- rate frequency with which lut is processed in Hz
- durationtotal length of stimulus flash in milliseconds. There is no guarantee that `duration %% length(lut)/rate == 0`. That is, the onus is on the user to ensure the duration is a multiple of the period of the stimuli.
- responseWindow maximum time ( $\geq 0$ ) in milliseconds to wait for a response from the onset of the stimulus presentation
- ... machine specific parameters

**Octopus 900**

x and y are in degrees, with precision to one decimal place recognised.

image is not possible on an Octopus 900.

lut is not possible on an Octopus 900. Stimulus is at 0 dB.

rate is in Hz, with precision to one decimal place recognised.

color is ignored. Use `opiSetBackground()` to alter stimulus color.

**Kowa AP7000**

Not supported.

**Compass**

Not implemented.

**See Also**

[opiSetBackground](#), [opiStaticStimulus](#), [opiKineticStimulus](#)

**Examples**

```
# A Size III flickering with a 10Hz square wave at location (7,7) with luminance
# 10 dB (HFA)
stim <- list(x=7, y=7, size=0.43, color="white",
            rate=20,          # one lut step per 50 ms
            lut=c(0,318),    # so one full lut per 100 ms == 10Hz
            duration=400,    # and 4 cycles per stimulus
            responseWindow=1500)
class(stim) <- "opiTemporalStimulus"
```

---

pixTodeg

*Convert pixels to degrees for machine 'machine'*

---

**Description**

Convert pixels to degrees for machine 'machine'

**Usage**

```
pixTodeg(xy, machine = "compass")
```

**Arguments**

xy                    a 2 element vector c(x,y) where x and y are in pixels  
 machine              "compass" or ...?

**Value**

xy converted to degrees of visual field with the usual conventions or NA if machine is unknown

**Examples**

```
pixTodeg(c(1000, 200), machine="compass") # c(1.290323, 24.516129) degrees
pixTodeg(c(1920/2, 1920/2)) # c(0,0) degrees
```

---

 QUESTP

*QUEST+*


---

### Description

An implementation of the Bayesian test procedure QUEST+ by AB Watson. This is mostly a translation of the MATLAB implementation by P Jones (see References). Its use is similar to ZEST. The objective is to estimate parameters of a function that defines the probability of responding stimuli. The steps are optimized based on entropy rather than the mean or the mode of the current pdfs.

### Usage

```

QUESTP(
    Fun,
    stimDomain,
    paramDomain,
    likelihoods = NULL,
    priors = NULL,
    stopType = "H",
    stopValue = 4,
    maxSeenLimit = 2,
    minNotSeenLimit = 2,
    minPresentations = 1,
    maxPresentations = 100,
    minInterStimInterval = NA,
    maxInterStimInterval = NA,
    verbose = 0,
    makeStim,
    ...
)

QUESTP.Prior(state, priors = NULL)

QUESTP.Likelihood(state)

QUESTP.start(
    Fun,
    stimDomain,
    paramDomain,
    likelihoods = NULL,
    priors = NULL,
    stopType = "H",
    stopValue = 4,
    maxSeenLimit = 2,
    minNotSeenLimit = 2,
    minPresentations = 1,

```



```

    maxPresentations = 100,
    makeStim,
    ...
)

getTargetStim(state)

QUESTP.step(state, nextStim = NULL)

QUESTP.stop(state)

QUESTP.final(state, Choice = "mean")

QUESTP.stdev(state, WhichP = NULL)

QUESTP.entropy(state)

```

### Arguments

Fun	Function to be evaluated, of the form <code>pseen = function(stim, param){...}</code> . Outputs a probability of seen.
stimDomain	Domain of values for the stimulus. Can be multi-dimensional (list, one element per dimension)
paramDomain	Domain of values for pdfs of the parameters in Fun. Can be multi-parametric (list, one element per parameter).
likelihoods	Pre-computed likelihoods if available (for QUESTP.start)
priors	Starting probability distributions for the parameter domains (list, one element per parameter)
stopType	N, for number of presentations; S, for standard deviation of the pdf; and H, for the entropy of the pdf (default).
stopValue	Value for number of presentations (stopType=N), standard deviation (stopType=S) or Entropy (stopType=H).
maxSeenLimit	Will terminate if maxStimulus value is seen this many times.
minNotSeenLimit	Will terminate if minStimulus value is not seen this many times.
minPresentations	Minimum number of presentations
maxPresentations	Maximum number of presentations regardless of stopType.
minInterStimInterval	If both minInterStimInterval and maxInterStimInterval are not NA, then between each stimuli there is a random wait period drawn uniformly between minInterStimInterval and maxInterStimInterval.
maxInterStimInterval	minInterStimInterval.

verbose	verbose=0 does nothing, verbose=1 stores pdfs for returning, and verbose=2 stores pdfs and also prints each presentaion.
makeStim	A function that takes a dB value and numPresentations and returns an OPI datatype ready for passing to opiPresent. See examples.
...	Extra parameters to pass to the opiPresent function
state	Current state of the QUESTP returned by QUESTP.start and QUESTP.step.
nextStim	A valid object for opiPresent to use as its nextStim.
Choice	How to compute final values in QUESTP.final ("mean","mode","median")
WhichP	Which parameter (numeric index) to monitor when calling QUESTP.stdev directly (returns max(stdev) if unspecified)

## Details

An implementation of the Bayesian test procedure QUEST+ by AB Watson. This is mostly a translation of the MATLAB implementation by P Jones (see References). Its use is similar to ZEST. The objective is to estimate parameters of a function that defines the probability of responding to stimuli. The steps are optimized based on entropy rather than the mean or the mode of the current pdfs.

The stimulus, parameter and response domain are separate and can be multidimensional. Each parameter has its own pdf. For evaluation, the pdfs are chained as a long vector (so no co-variances are considered). More complex functions will require larger combined pdfs and longer computation time for likelihoods and updating at each step. In these cases, it is recommended to pre-calculate the likelihoods using QUESTP.Likelihood and store them.

The function to be fitted needs to output a probability of seen (i.e. `pseen = function(stim, param){...}`) and must take `stim` and `param` as inputs. `stim` is a vector with length = number of stimulus dimensions (in simple one-dimensional cases, the intensity in dB). `param` is a vector with length = number of parameters to be fitted in `Fun`.

For example, QUEST+ can fit a Gaussian psychometric function with `stimDomain = {0, 1, ..., 39, 40}` dB and `paramDomain = ({0, 1, ..., 39, 40}; {0.5, 1, ..., 5.5, 6})` dB for the mean and standard deviation respectively. A standard ZEST procedure can be replicated by setting `stimDomain = {0, 1, ..., 39, 40}` dB and `paramDomain = ({0, 1, ..., 39, 40}; {1})` dB, i.e. by setting the `stimDomain = paramDomain` for the mean and by having a static standard deviation = 1 dB. Note however that the stimulus selection is always based on entropy and not on the mean/mode of the current pdf. See examples below

Note this function will repeatedly call `opiPresent` for a stimulus until `opiPresent` returns NULL (ie no error occurred).

The `checkFixationOK` function is called (if present in `stim` made from `makeStim`) after each presentation, and if it returns FALSE, the pdf for that location is not changed (ie the presentation is ignored), but the `stim`, number of presentations etc is recorded in the `state`.

If more than one QUESTP is to be interleaved (for example, testing multiple locations), then the `QUESTP.start`, `QUESTP.step`, `QUESTP.stop` and `QUESTP.final` calls can maintain the state of the QUESTP after each presentation, and should be used. If only a single QUESTP is required, then the simpler QUESTP can be used, which is a wrapper for the four functions that maintain state. See examples below.

**Value****Single location:**

QUESTP returns a list containing

- `npres` Total number of presentations used.
- `respSeq` Response sequence stored as a data frame: column 1 is a string identified of a (potentially) multidimensional stimulus values of stimuli (dimensions chained into a string), column 2 is 1/0 for seen/not-seen, column 3 is fixated 1/0 (always 1 if `checkFixationOK` not present in stim objects returned from `makeStim`). All additional columns report each stimulus dimension, one for each row.
- `pdfs` If `verbose` is bigger than 0, then this is a list of the pdfs used for each presentation, otherwise NULL.
- `final` The mean (default, strongly suggested)/median/mode of the parameters' pdf, depending on `Choice`.
- `opiResp` A list of responses received from each successful call to `opiPresent` within QUESTP.

**Multiple locations:**

QUESTP.`start` returns a list that can be passed to QUESTP.`step`, QUESTP.`stop`, and QUESTP.`final`. It represents the state of a QUESTP at a single location at a point in time and contains the following.

- `name` QUESTP
- A copy of all of the parameters supplied to QUESTP.`start`: `stimDomain`, `paramDomain`, `likelihoods`, `priors`, `stopType`, `stopValue`, `maxSeenLimit`, `minNotSeenLimit`, `minPresentations`, `maxPresentations`, `makeStim`, and `opiParams`.
- `pdf` Current pdf: vector of probabilities, collating all parameter domains.
- `priorsP` List of starting pdfs, one for each parameter.
- `numPresentations` The number of times QUESTP.`step` has been called on this state.
- `stimuli` A vector containing the stimuli used at each call of QUESTP.`step`.
- `responses` A vector containing the responses received at each call of QUESTP.`step`.
- `responseTimes` A vector containing the response times received at each call of QUESTP.`step`.
- `fixated` A vector containing TRUE/FALSE if fixation was OK according to `checkFixationOK` for each call of QUESTP.`step` (defaults to TRUE if `checkFixationOK` not present).
- `opiResp` A list of responses received from each call to `opiPresent` within QUESTP.`step`.

QUESTP.`step` returns a list containing

- `state` The new state after presenting a stimuli and getting a response.
- `resp` The return from the `opiPresent` call that was made.

QUESTP.`stop` returns TRUE if the QUESTP has reached its stopping criteria, and FALSE otherwise.

QUESTP.`final` returns an estimate of parameters based on state. If `state$Choice` is `mean` then the mean is returned (the only one that really makes sense for QUEST+). If `state$Choice` is `mode` then the mode is returned. If `state$Choice` is `median` then the median is returned.

**References**

Andrew B. Watson; QUEST+: A general multidimensional Bayesian adaptive psychometric method. *Journal of Vision* 2017;17(3):10. doi: <https://doi.org/10.1167/17.3.10>.

Jones, P. R. (2018). QuestPlus: a MATLAB implementation of the QUEST+ adaptive psychometric method, *Journal of Open Research Software*, 6(1):27. doi: <http://doi.org/10.5334/jors.195>

A. Turpin, P.H. Artes and A.M. McKendrick "The Open Perimetry Interface: An enabling tool for clinical visual psychophysics", *Journal of Vision* 12(11) 2012.

### See Also

[dbTocd](#), [opiPresent](#)

### Examples

```
chooseOpi("SimHenson")
if(!is.null(opiInitialize(type="C", cap=6)))
  stop("opiInitialize failed")

#####
# This section is for single location QUESTP
# This example fits a FoS curve
# Note: only fitting threshold and slope,
# modify the domain for FPR and FNR to fit those as well
#####
# Stimulus is Size III white-on-white as in the HFA
makeStim <- function(db, n) {
  s <- list(x=9, y=9, level=dbTocd(db), size=0.43, color="white",
           duration=200, responseWindow=1500, checkFixationOK=NULL)
  class(s) <- "opiStaticStimulus"
  return(s)
}

#True parameters (variability is determined according to Henson et al. based on threshold)
loc <- list(threshold = 20, fpr = 0.05, fnr = 0.05)

#Function to fit (Gaussian psychometric function)
pSeen <- function(x, params){return(params[3] +
                                   (1 - params[3] - params[4]) *
                                   (1 - pnorm(x, params[1], params[2])))}

#QUEST+
QP <- QUESTP(Fun = pSeen,
             stimDomain = list(0:50),
             paramDomain = list(seq(0, 40, 1), #Domain for the 50% threshold (Mean)
                               seq(.5, 8, .5), #Domain for the slope (SD)
                               seq(0.05, 0.05, 0.05), #Domain for the FPR (static)
                               seq(0.05, 0.05, 0.05)), #Domain for the FNR (static)
             stopType="H", stopValue=4, maxPresentations=500,
             makeStim = makeStim,
             tt=loc$threshold, fpr=loc$fpr, fnr=loc$fnr,
             verbose = 2)

#Plots results
#Henson's FoS function (as implemented in OPI - ground truth)
HensFunction <- function(Th){
  SD <- exp(-0.081*Th + 3.27)
```

```

    SD[SD > 6] <- 6
    return(SD)
}

#Stimulus domain
dB_Domain <- 0:50
FoS <- pSeen(dB_Domain, params = QP$final) # Estimated FoS
FoS_GT <- pSeen(dB_Domain, params = c(loc$threshold, HensFunction(loc$threshold),
                                     loc$fpr, loc$fnr)) #Ground truth FoS (based on Henson et al.)

#Plot (seen stimuli at the top, unseen stimuli at the bottom)
plot(dB_Domain, FoS_GT, type = "l", ylim = c(0, 1), xlab = "dB", ylab = "% seen", col = "blue")
lines(dB_Domain, FoS, col = "red")
points(QP$respSeq$stimuli, QP$respSeq$responses, pch = 16,
       col = rgb(red = 1, green = 0, blue = 0, alpha = 0.1))
legend("top", inset = c(0, -.2), legend = c("True", "Estimated", "Stimuli"),
      col=c("blue", "red", "red"), lty=c(1,1,0),
      pch = c(16, 16, 16), pt.cex = c(0, 0, 1),
      horiz = TRUE, xpd = TRUE, xjust = 0)

if (!is.null(opiClose()))
  warning("opiClose() failed")

chooseOpi("SimHenson")
if(!is.null(opiInitialize(type="C", cap=6)))
  stop("opiInitialize failed")

#####
# This section is for single location QUESTP
# This example shows that QUEST+ can replicate a ZEST procedure
# by fitting a FoS curve with fixed Slope, FPR and FNR
# Compared with ZEST
# Note that QUEST+ should be marginally more efficient in selecting
# the most informative stimulus
#####
# Stimulus is Size III white-on-white as in the HFA
makeStim <- function(db, n) {
  s <- list(x=9, y=9, level=dbTocd(db), size=0.43, color="white",
           duration=200, responseWindow=1500, checkFixationOK=NULL)
  class(s) <- "opiStaticStimulus"
  return(s)
}

#True parameters (variability is determined according to Henson et al. based on threshold)
loc <- list(threshold = 30, fpr = 0.03, fnr = 0.03)

#Function to fit (Gaussian psychometric function - Fixed slope (same as default likelihood in ZEST))
pSeen <- function(domain, tt){0.03+(1-0.03-0.03)*(1-pnorm(domain,tt,1))}

# ZEST-like QUEST+ procedure
QP <- QUESTP(Fun = pSeen,
```

```

stimDomain = list(0:40),
paramDomain = list(seq(0, 40, 1)),
stopType="S", stopValue=1.5, maxPresentations=500,
makeStim = makeStim,
tt=loc$threshold, fpr=loc$fpr, fnr=loc$fnr,
verbose = 2)

# ZEST
ZE <- ZEST(domain = 0:40,
stopType="S", stopValue=1.5, maxPresentations=500,
makeStim = makeStim,
tt=loc$threshold, fpr=loc$fpr, fnr=loc$fnr,
verbose = 2)

#Plots results
#Henson's FoS function (as implemented in OPI - ground truth)
HensFunction <- function(Th){
SD <- exp(-0.081*Th + 3.27)
SD[SD > 6] <- 6
return(SD)
}

#Stimulus domain
dB_Domain <- 0:50
FoS_QP <- pSeen(domain = dB_Domain, tt = QP$final) # Estimated FoS
FoS_ZE <- pSeen(domain = dB_Domain, tt = ZE$final) # Estimated FoS

#Plot (seen stimuli at the top, unseen stimuli at the bottom)
plot(dB_Domain, FoS_QP, type = "l", ylim = c(0, 1), xlab = "dB", ylab = "% seen", col = "blue")
lines(dB_Domain, FoS_ZE, col = "red")
points(QP$respSeq$stimuli, QP$respSeq$responses, pch = 16,
col = rgb(red = 0, green = 0, blue = 1, alpha = 0.5))
points(ZE$respSeq[1,], ZE$respSeq[2,], pch = 16,
col = rgb(red = 1, green = 0, blue = 0, alpha = 0.5))
legend("bottomleft", legend = c("QUEST+", "ZEST", "Stimuli QUEST+", "Stimuli ZEST"),
col=c("blue", "red", "blue", "red"), lty=c(1,1,0,0),
pch = c(16, 16, 16, 16), pt.cex = c(0, 0, 1, 1),
horiz = FALSE, xpd = TRUE, xjust = 0)
abline(v = loc$threshold, lty = "dashed")

if (!is.null(opiClose()))
warning("opiClose() failed")

chooseOpi("SimHenson")
if(!is.null(opiInitialize(type="C", cap=6)))
stop("opiInitialize failed")

#####
# This section is for single location QUESTP
# This example fits a broken stick spatial summation function
# with a multi-dimensional stimulus (varying in size and intensity).
# Stimulus sizes are limited to GI, GII, GIII, GIV and GV.

```

```

# The example also shows how to use a helper function to
# simulate responses to multi-dimensional stimuli
# (here, the simulated threshold varies based on stimulus size)
#####
makeStim <- function(stim, n) {
  s <- list(x=9, y=9, level=dbTocd(stim[1]), size=stim[2], color="white",
           duration=200, responseWindow=1500, checkFixationOK=NULL)
  class(s) <- "opiStaticStimulus"
  return(s)
}

# Helper function for true threshold (depends on log10(stimulus size),
# diameter assumed to be the second element of stim vector)
ttHelper_SS <- function(location) { # returns a function of (stim)
  ff <- function(stim) stim

  body(ff) <- substitute(
    {return(SensF(log10(pi*(stim[2]/2)^2), c(location$Int1, location$Int2, location$Slo2))})
  )
  return(ff)
}

# Function of sensitivity vs SSize (log10(stimulus area))
SensF <- function(SSize, params){
  Sens <- numeric(length(SSize))
  for (i in 1:length(SSize)){
    Sens[i] <- min(c(params[1] + 10*SSize[i], params[2] + params[3]*SSize[i]))
  }
  Sens[Sens < 0] <- 0
  return(Sens)
}

Sizes <- c(0.1, 0.21, 0.43, 0.86, 1.72)

#True parameters (variability is determined according to Henson et al. based on threshold)
loc <- list(Int1 = 32, Int2 = 28, Slo2 = 2.5, fpr = 0.05, fnr = 0.05, x = 9, y = 9)

# Function to fit (probability of seen given a certain stimulus intensity and size,
# for different parameters)
pSeen <- function(stim, params){

  Th <- SensF(log10(pi*(stim[2]/2)^2), params)

  return(0.03 +
         (1 - 0.03 - 0.03) *
         (1 - pnorm(stim[1], Th, 1)))
}

set.seed(111)
#QUEST+ - takes some time to calculate likelihoods
## Not run:

```

```

QP <- QUESTP(Fun = pSeen,
             stimDomain = list(0:50, Sizes),
             paramDomain = list(seq(0, 40, 1), # Domain for total summation intercept
                               seq(0, 40, 1), # Domain for partial summation intercept
                               seq(0, 3, 1)), # Domain for partial summation slope
             stopType="H", stopValue=1, maxPresentations=500,
             makeStim = makeStim,
             ttHelper=ttHelper_SS(loc), tt = 30,
             fpr=loc$fpr, fnr=loc$fnr,
             verbose = 2)

#Stimulus sizes
G <- log10(c(pi*(0.1/2)^2, pi*(0.21/2)^2, pi*(0.43/2)^2, pi*(0.86/2)^2, pi*(1.72/2)^2));
SizesP <- seq(min(G), max(G), .05)

# True and estimated response
Estim_Summation <- SensF(SizesP, params = QP$final) # Estimated spatial summation
GT_Summation <- SensF(SizesP, params = c(loc$Int1, loc$Int2, loc$Slo2)) # True spatial summation

#Plot
plot(10^SizesP, GT_Summation, type = "l", ylim = c(0, 40), log = "x",
     xlab = "Stimulus area (deg^2)", ylab = "Sensitivity (dB)", col = "blue")
lines(10^SizesP, Estim_Summation, col = "red")
points(pi*(QP$respSeq$stimuli.2/2)^2, QP$respSeq$stimuli.1, pch = 16,
       col = rgb(red = 1, green = 0, blue = 0, alpha = 0.3))
legend("top", inset = c(0, -.2), legend = c("True", "Estimated", "Stimuli"),
      col=c("blue", "red", "red"), lty=c(1,1,0),
      pch = c(16, 16, 16), pt.cex = c(0, 0, 1),
      horiz = TRUE, xpd = TRUE, xjust = 0)

## End(Not run)

if (!is.null(opiClose()))
  warning("opiClose() failed")

```

---

RtDbUnits

*Response times to white-on-white Goldmann Size III targets for 12 subjects in dB units*

---

### Description

Response times to white-on-white Goldmann Size III targets for 12 subjects. The second column is the distance of the stimuli from measured threshold in HFA dB units. The threshold was determined by post-hoc fit of FoS curves to the data.

### Usage

RtDbUnits



**Format**

An object of class `data.frame` with 30620 rows and 3 columns.

**Details**

A data frame with 30620 observations on the following 3 variables.

- Rt Reaction time in ms.
- Dist Distance of stimuli from threshold in dB.
- Person Identifier of each subject.

**References**

A.M. McKendrick, J. Denniss and A. Turpin. "Response times across the visual field: empirical observations and application to threshold determination". *Vision Research* 101 2014.

**See Also**

[RtSigmaUnits](#)

---

RtSigmaUnits

*Response times to white-on-white Goldmann Size III targets for 12 subjects in sigma units*

---

**Description**

Response times to white-on-white Goldmann Size III targets for 12 subjects. The second column is the distance of the stimuli from measured threshold in 'sigma' units. The threshold was determined by post-hoc fit of a cumulative gaussian FoS curve to the data for each location and subject. Sigma is the standard deviation of the fitted FoS.

**Usage**

```
RtSigmaUnits
```

**Format**

An object of class `data.frame` with 30620 rows and 3 columns.

**Details**

A data frame with 30620 observations on the following 3 variables.

- Rt Reaction time in ms.
- Dist Distance of stimuli from threshold in sigma units.
- Person Identifier of each subject.

## References

A.M. McKendrick, J. Denniss and A. Turpin. "Response times across the visual field: empirical observations and application to threshold determination". *Vision Research* 101 2014.

## See Also

[RtDbUnits](#)

---

ZEST

*ZEST*

---

## Description

An implementation of the Bayesian test procedures of King-Smith et al. and Watson and Pelli. Note that we use the term pdf throughout as in the original paper, even though they are discrete probability functions in this implementation.

## Usage

```
ZEST(
  domain = 0:40,
  prior = rep(1/length(domain), length(domain)),
  likelihood = sapply(domain, function(tt) 0.03 + (1 - 0.03 - 0.03) * (1 -
    stats::pnorm(domain, tt, 1))),
  stopType = "S",
  stopValue = 1.5,
  minStimulus = utils::head(domain, 1),
  maxStimulus = utils::tail(domain, 1),
  maxSeenLimit = 2,
  minNotSeenLimit = 2,
  maxPresentations = 100,
  minInterStimInterval = NA,
  maxInterStimInterval = NA,
  verbose = 0,
  makeStim,
  stimChoice = "mean",
  ...
)
```

```
ZEST.start(
  domain = 0:40,
  prior = rep(1/length(domain), length(domain)),
  likelihood = sapply(domain, function(tt) 0.03 + (1 - 0.03 - 0.03) * (1 -
    stats::pnorm(domain, tt, 1))),
  stopType = "S",
  stopValue = 1.5,
  minStimulus = utils::head(domain, 1),
```

```

    maxStimulus = utils::tail(domain, 1),
    maxSeenLimit = 2,
    minNotSeenLimit = 2,
    maxPresentations = 100,
    makeStim,
    stimChoice = "mean",
    ...
)

ZEST.step(state, nextStim = NULL, fixedStimValue = NA)

ZEST.stop(state)

ZEST.final(state)

```

### Arguments

domain	Vector of values over which pdf is kept.
prior	Starting probability distribution over domain. Same length as domain.
likelihood	Matrix where $\text{likelihood}[s, t]$ is likelihood of seeing $s$ given $t$ is the true threshold. That is, $\text{Pr}(s t)$ where $s$ and $t$ are indexes into domain.
stopType	N, for number of presentations; S, for standard deviation of the pdf; and H, for the entropy of the pdf.
stopValue	Value for number of presentations ( $\text{stopType}=\text{N}$ ), standard deviation ( $\text{stopType}=\text{S}$ ) or Entropy ( $\text{stopType}=\text{H}$ ).
minStimulus	The smallest stimuli that will be presented. Could be different from $\text{domain}[1]$ .
maxStimulus	The largest stimuli that will be presented. Could be different from $\text{tail}(\text{domain}, 1)$ .
maxSeenLimit	Will terminate if maxStimulus value is seen this many times.
minNotSeenLimit	Will terminate if minStimulus value is not seen this many times.
maxPresentations	Maximum number of presentations regardless of stopType.
minInterStimInterval	If both minInterStimInterval and maxInterStimInterval are not NA, then between each stimuli there is a random wait period drawn uniformly between minInterStimInterval and maxInterStimInterval.
maxInterStimInterval	minInterStimInterval.
verbose	verbose=0 does nothing, verbose=1 stores pdfs for returning, and verbose=2 stores pdfs and also prints each presentation.
makeStim	A function that takes a dB value and numPresentations and returns an OPI datatype ready for passing to opiPresent. See examples.
stimChoice	A true ZEST procedure uses the "mean" of the current pdf as the stimulus, but "median" and "mode" (as used in a QUEST procedure) are provided for your enjoyment.

...	Extra parameters to pass to the <code>opiPresent</code> function
<code>state</code>	Current state of the ZEST returned by <code>ZEST.start</code> and <code>ZEST.step</code> .
<code>nextStim</code>	A valid object for <code>opiPresent</code> to use as its <code>nextStim</code> .
<code>fixedStimValue</code>	A number in <code>state\$domain</code> that, is <code>!is.na</code> , will be used as the stimulus value overriding <code>state\$minStimulus</code> , <code>state\$maxStimulus</code> and <code>state\$stimChoice</code> .

## Details

This is an implementation of King-Smith et al.'s ZEST procedure and Watson and Pelli's QUEST procedure. All presentations are rounded to an element of the supplied domain.

Note this function will repeatedly call `opiPresent` for a stimulus until `opiPresent` returns NULL (ie no error occurred).

The `checkFixationOK` function is called (if present in `stim` made from `makeStim`) after each presentation, and if it returns FALSE, the pdf for that location is not changed (ie the presentation is ignored), but the `stim`, number of presentations etc is recorded in the state.

If more than one ZEST is to be interleaved (for example, testing multiple locations), then the `ZEST.start`, `ZEST.step`, `ZEST.stop` and `ZEST.final` calls can maintain the state of the ZEST after each presentation, and should be used. If only a single ZEST is required, then the simpler ZEST can be used, which is a wrapper for the four functions that maintain state. See examples below.

## Value

### Single location:

ZEST returns a list containing

- `npres` Total number of presentations used.
- `respSeq` Response sequence stored as a matrix: row 1 is dB values of stimuli, row 2
- `is` 1/0 for seen/not-seen, row 3 is fixated 1/0 (always 1 if `checkFixationOK` not
- present in `stim` objects returned from `makeStim`).
- `pdfs` If `verbose` is bigger than 0, then this is a list of the pdfs used for each presentation, otherwise NULL.
- `final` The mean/median/mode of the final pdf, depending on `stimChoice`, which is the determined threshold.
- `opiResp` A list of responses received from each successful call to `opiPresent` within ZEST.

### Multiple locations:

`ZEST.start` returns a list that can be passed to `ZEST.step`, `ZEST.stop`, and `ZEST.final`. It represents the state of a ZEST at a single location at a point in time and contains the following.

- `name` ZEST.
- `pdf` Current pdf: vector of probabilities the same length as domain.
- `numPresentations` The number of times `ZEST.step` has been called on this state.
- `stimuli` A vector containing the stimuli used at each call of `ZEST.step`.
- `responses` A vector containing the responses received at each call of `ZEST.step`.
- `responseTimes` A vector containing the response times received at each call of `ZEST.step`.

- `fixated` A vector containing TRUE/FALSE if fixation was OK according to `checkFixationOK` for each call of `ZEST.step` (defaults to TRUE if `checkFixationOK` not present).
- `opiResp` A list of responses received from each call to `opiPresent` within `ZEST.step`.
- A copy of all of the parameters supplied to `ZEST.start`: `domain`, `likelihood`, `stopType`, `stopValue`, `minStimulus`, `maxStimulus`, `maxSeenLimit`, `minNotSeenLimit`, `maxPresentations`, `makeStim`, `stimChoice`, `currSeenLimit`, `currNotSeenLimit`, and `opiParams`.

`ZEST.step` returns a list containing \* `state` The new state after presenting a stimuli and getting a response. \* `resp` The return from the `opiPresent` call that was made.

`ZEST.stop` returns TRUE if the ZEST has reached its stopping criteria, and FALSE otherwise.

`ZEST.final` returns an estimate of threshold based on state. If `state$stimChoice` is mean then the mean is returned. If `state$stimChoice` is mode then the mode is returned. If `state$stimChoice` is median then the median is returned.

## References

P.E. King-Smith, S.S. Grigsby, A.J. Vingrys, S.C. Benes, and A. Supowit. "Efficient and Unbiased Modifications of the QUEST Threshold Method: Theory, Simulations, Experimental Evaluation and Practical Implementation", *Vision Research* 34(7) 1994. Pages 885-912.

A.B. Watson and D.G. Pelli. "QUEST: A Bayesian adaptive psychophysical method", *Perception and Psychophysics* 33 1983. Pages 113-120.

A. Turpin, P.H. Artes and A.M. McKendrick "The Open Perimetry Interface: An enabling tool for clinical visual psychophysics", *Journal of Vision* 12(11) 2012.

## See Also

[dbToccd](#), [opiPresent](#)

## Examples

```
chooseOpi("SimHenson")
if(!is.null(opiInitialize(type="C", cap=6)))
  stop("opiInitialize failed")

#####
# This section is for single location ZESTs
#####
# Stimulus is Size III white-on-white as in the HFA
makeStim <- function(db, n) {
  s <- list(x=9, y=9, level=dbToccd(db), size=0.43, color="white",
           duration=200, responseWindow=1500, checkFixationOK=NULL)
  class(s) <- "opiStaticStimulus"
  return(s)
}

repp <- function(...) sapply(1:50, function(i) ZEST(makeStim=makeStim, ...))
a <- repp(stopType="H", stopValue= 3, verbose=0, tt=30, fpr=0.03)
b <- repp(stopType="S", stopValue=1.5, verbose=0, tt=30, fpr=0.03)
c <- repp(stopType="S", stopValue=2.0, verbose=0, tt=30, fpr=0.03)
d <- repp(stopType="N", stopValue= 50, verbose=0, tt=30, fpr=0.03)
```

```

e <- repp(prior=dnorm(0:40,m=0,s=5), tt=30, fpr=0.03)
f <- repp(prior=dnorm(0:40,m=10,s=5), tt=30, fpr=0.03)
g <- repp(prior=dnorm(0:40,m=20,s=5), tt=30, fpr=0.03)
h <- repp(prior=dnorm(0:40,m=30,s=5), tt=30, fpr=0.03)

layout(matrix(1:2,1,2))
boxplot(lapply(list(a,b,c,d,e,f,g,h), function(x) unlist(x["final",])))
boxplot(lapply(list(a,b,c,d,e,f,g,h), function(x) unlist(x["npres",])))

#####
# This section is for multiple ZESTs
#####
makeStimHelper <- function(db,n, x, y) { # returns a function of (db,n)
  ff <- function(db, n) db+n
  body(ff) <- substitute({
    s <- list(x=x, y=y, level=dbTocd(db), size=0.43, color="white",
              duration=200, responseWindow=1500, checkFixationOK=NULL)
    class(s) <- "opiStaticStimulus"
    return(s)
  }, list(x=x,y=y))
  return(ff)
}

# List of (x, y, true threshold) triples
locations <- list(c(9,9,30), c(-9,-9,32), c(9,-9,31), c(-9,9,33))

# Setup starting states for each location
states <- lapply(locations, function(loc) {
  ZEST.start(
    domain=-5:45,
    minStimulus=0,
    maxStimulus=40,
    makeStim=makeStimHelper(db,n,loc[1],loc[2]),
    stopType="S", stopValue= 1.5, tt=loc[3], fpr=0.03, fnr=0.01)})

# Loop through until all states are "stop"
while(!all(st <- unlist(lapply(states, ZEST.stop)))) {
  i <- which(!st) # choose a random,
  i <- i[runif(1, min=1, max=length(i))] # unstopped state
  r <- ZEST.step(states[[i]]) # step it
  states[[i]] <- r$state # update the states
}

finals <- lapply(states, ZEST.final) # get final estimates of threshold
for(i in 1:length(locations)) {
  #cat(sprintf("Location (%+2d,%+2d) ",locations[[i]][1], locations[[i]][2]))
  #cat(sprintf("has threshold %4.2f\n", finals[[i]]))
}

if (!is.null(opiClose()))
  warning("opiClose() failed")

```

# Index

- \* **datasets**
  - .OpiEnv, 2
  - opi.implementations, 16
- \* **dataset**
  - RtDbUnits, 56
  - RtSigmaUnits, 57
  - .OpiEnv, 2
- cdTodb, 3
- chooseOPI (chooseOpi), 3
- chooseOpi, 3, 16, 18, 25, 34, 37, 43
- compass.opiClose (opiClose), 16
- compass.opiInitialize (opiInitialize), 19
- compass.opiPresent (opiPresent), 28
- compass.opiQueryDevice (opiQueryDevice), 35
- compass.opiSetBackground (opiSetBackground), 37
- daydream.opiClose (opiClose), 16
- daydream.opiInitialize (opiInitialize), 19
- daydream.opiPresent (opiPresent), 28
- daydream.opiQueryDevice (opiQueryDevice), 35
- daydream.opiSetBackground (opiSetBackground), 37
- dbTocd, 5, 7, 10, 14, 52, 61
- degTopix, 5
- display.opiClose (opiClose), 16
- display.opiInitialize (opiInitialize), 19
- display.opiPresent (opiPresent), 28
- display.opiQueryDevice (opiQueryDevice), 35
- display.opiSetBackground (opiSetBackground), 37
- fourTwo.final (fourTwo.start), 6
- fourTwo.start, 6, 10
- fourTwo.step (fourTwo.start), 6
- fourTwo.stop (fourTwo.start), 6
- FT, 7, 8
- getTargetStim (QUESTP), 48
- imo.opiClose (opiClose), 16
- imo.opiInitialize (opiInitialize), 19
- imo.opiPresent (opiPresent), 28
- imo.opiQueryDevice (opiQueryDevice), 35
- imo.opiSetBackground (opiSetBackground), 37
- kowaAP7000.opiClose (opiClose), 16
- kowaAP7000.opiInitialize (opiInitialize), 19
- kowaAP7000.opiPresent (opiPresent), 28
- kowaAP7000.opiQueryDevice (opiQueryDevice), 35
- kowaAP7000.opiSetBackground (opiSetBackground), 37
- MOCS, 11
- octo600.opiClose (opiClose), 16
- octo600.opiInitialize (opiInitialize), 19
- octo600.opiPresent (opiPresent), 28
- octo600.opiQueryDevice (opiQueryDevice), 35
- octo600.opiSetBackground (opiSetBackground), 37
- octo900.opiClose (opiClose), 16
- octo900.opiInitialize (opiInitialize), 19
- octo900.opiPresentF310 (opiPresent), 28
- octo900.opiQueryDevice (opiQueryDevice), 35
- octo900.opiSetBackground (opiSetBackground), 37

opi.implementations, 16  
 opiClose, 16, 25  
 opiDistributor (opi.implementations), 16  
 opiGetParams, 18  
 opiInitialise, 30  
 opiInitialise (opiInitialize), 19  
 opiInitialize, 19, 34  
 opiKineticStimulus, 26, 29, 34, 45, 47  
 opiPresent, 7, 10, 14, 25, 28, 29, 30, 34, 52, 61  
 opiQueryDevice, 35  
 opiSetBackground, 25, 28, 37, 45, 47  
 opiStaticStimulus, 28, 29, 34, 44, 47  
 opiTemporalStimulus, 28, 29, 34, 45, 46

PhoneHMD.opiClose (opiClose), 16  
 PhoneHMD.opiInitialize (opiInitialize), 19  
 PhoneHMD.opiPresent (opiPresent), 28  
 PhoneHMD.opiQueryDevice (opiQueryDevice), 35  
 PhoneHMD.opiSetBackground (opiSetBackground), 37  
 pixTodeg, 47

QUESTP, 48

RtDbUnits, 56, 58  
 RtSigmaUnits, 57, 57

simG.opiClose (opiClose), 16  
 simG.opiInitialize (opiInitialize), 19  
 simG.opiPresent (opiPresent), 28  
 simG.opiQueryDevice (opiQueryDevice), 35  
 simG.opiSetBackground (opiSetBackground), 37  
 simH.opiInitialize (opiInitialize), 19  
 simH.opiPresent (opiPresent), 28  
 simH.opiSetBackground (opiSetBackground), 37  
 simH\_RT.opiInitialize (opiInitialize), 19  
 simH\_RT.opiPresent (opiPresent), 28  
 simH\_RT.opiSetBackground (opiSetBackground), 37  
 simNo.opiPresent (opiPresent), 28  
 simNo.opiSetBackground (opiSetBackground), 37  
 simYes.opiPresent (opiPresent), 28

simYes.opiSetBackground (opiSetBackground), 37

ZEST, 58