

# Package ‘ieegio’

September 16, 2024

**Title** File IO for Intracranial Electroencephalography

**Version** 0.0.1

**Language** en-US

**Encoding** UTF-8

**Description** Integrated toolbox supporting common file formats used for intracranial  
Electroencephalography (iEEG) and deep-brain stimulation (DBS) study.

**URL** <http://dipterix.org/ieegio/>

**BugReports** <https://github.com/dipterix/ieegio/issues>

**License** MIT + file LICENSE

**RoxygenNote** 7.3.2

**Imports** data.table (>= 1.16.0), digest, fastmap, filearray (>= 0.1.8),  
freesurferformats, fs, fst (>= 0.9.0), gifti (>= 0.8.0),  
grDevices, hdf5r, jsonlite, oro.nifti, R.matlab (>= 3.7.0), R6,  
readNSx (>= 0.0.5), rpyANTs (>= 0.0.3), stringr, utils, yaml

**Suggests** reticulate, rgl, rpymat (>= 0.1.7), RNifti (>= 1.7.0), xml2,  
knitr, rmarkdown, tools

**VignetteBuilder** knitr

**NeedsCompilation** no

**Author** Zhengjia Wang [aut, cre] (<<https://orcid.org/0000-0001-5629-1116>>)

**Maintainer** Zhengjia Wang <dipterix.wang@gmail.com>

**Repository** CRAN

**Date/Publication** 2024-09-16 14:20:02 UTC

## Contents

ieegio_sample_data . . . . .	2
imaging-surface . . . . .	3
imaging-volume . . . . .	5
io_h5_valid . . . . .	7
io_read_h5 . . . . .	9

io_write_h5 . . . . .	10
LazyH5 . . . . .	11
low-level-read-write . . . . .	14
NWBHDF5IO . . . . .	17
plot.ieegio_surface . . . . .	21
plot.ieegio_volume . . . . .	23
pynwb_module . . . . .	25
read_bci2000 . . . . .	26
read_brainvis . . . . .	27
read_edf . . . . .	28
read_nsx . . . . .	29
read_nwb . . . . .	30
SignalDataCache . . . . .	33

<b>Index</b>	<b>35</b>
--------------	-----------

---

ieegio\_sample\_data     *Download sample files*

---

## Description

Download sample files

## Usage

```
ieegio_sample_data(file, test = FALSE, cache_ok = TRUE)
```

## Arguments

file	file to download; set to NULL to view all possible files
test	test whether the sample file exists instead of downloading them; default is FALSE
cache_ok	whether to use cache

## Value

When test is false, returns downloaded file path (character); when test is true, returns whether the expected sample exists (logical).

## Examples

```
# list available files
ieegio_sample_data()

# check if file edfPlusD.edf exists
ieegio_sample_data("edfPlusD.edf", test = TRUE)

## Not run:
ieegio_sample_data("edfPlusD.edf")
```

```
## End(Not run)
```

---

imaging-surface	<i>Read and write surface files</i>
-----------------	-------------------------------------

---

## Description

Supports surface geometry, annotation, measurement, and time-series data. Please use the high-level function `read_surface`, which calls other low-level functions internally.

## Usage

```
read_surface(file, format = "auto", type = NULL, ...)

write_surface(
  x,
  con,
  format = c("gifti", "freesurfer"),
  type = c("geometry", "annotations", "measurements", "color", "time_series"),
  ...
)

io_read_fs(
  file,
  type = c("geometry", "annotations", "measurements"),
  format = "auto",
  name = basename(file),
  ...
)

io_read_gii(file)

io_write_gii(x, con, ...)
```

## Arguments

file, con	path the file
format	format of the file, see 'Arguments' section in <code>read.fs.surface</code> (when file type is 'geometry') and <code>read.fs.curv</code> (when file type is 'measurements')
type	type of the data; ignored if the file format is 'GIFTI'. For 'FreeSurfer' files, supported types are 'geometry' contains positions of mesh vertex nodes and face indices;

'annotations' annotation file (usually with file extension 'annot') containing a color look-up table and an array of color keys. These files are used to display discrete values on the surface such as brain atlas;  
 'measurements' measurement file such as 'sulc' and 'curv' files, containing numerical values (often with continuous domain) for each vertex node  
 ...  
 for `read_surface`, the arguments will be passed to `io_read_fs` if the file is a 'FreeSurfer' file.  
 x surface (geometry, annotation, measurement) data  
 name name of the data; default is the file name

## Value

A surface object container

## Examples

```

library(ieegio)

# geometry
geom_file <- "gifti/GzipBase64/sujet01_Lwhite.surf.gii"

# measurements
shape_file <- "gifti/GzipBase64/sujet01_Lwhite.shape.gii"

# time series
ts_file <- "gifti/GzipBase64/fmri_sujet01_Lwhite_projection.time.gii"

if(ieegio_sample_data(geom_file, test = TRUE)) {

  geometry <- read_surface(ieegio_sample_data(geom_file))
  print(geometry)

  measurement <- read_surface(ieegio_sample_data(shape_file))
  print(measurement)

  time_series <- read_surface(ieegio_sample_data(ts_file))
  print(time_series)

  # merge measurement & time_series into geometry
  merged <- merge(geometry, measurement, time_series)
  print(merged)

  # make sure you install `rgl` package
  plot(merged, name = c("measurements", "Shape001"))

  plot(merged, name = "time_series",
       slice_index = c(1, 11, 21, 31))

}

```

---

imaging-volume	<i>Read and write volume data</i>
----------------	-----------------------------------

---

## Description

Read and write volume data ('MRI', 'CT', etc.) in 'NIfTI' or 'MGH' formats. Please use `read_volume` and `write_volume` for high-level function. These functions will call other low-level functions internally.

## Usage

```
read_volume(file, header_only = FALSE, format = c("auto", "nifti", "mgh"), ...)

write_volume(x, con, format = c("auto", "nifti", "mgh"), ...)

io_read_mgz(file, header_only = FALSE)

io_write_mgz(x, con, ...)

## S3 method for class 'ieegio_volume'
io_write_mgz(x, con, ...)

## S3 method for class 'ieegio_mgh'
io_write_mgz(x, con, ...)

## S3 method for class 'nifti'
io_write_mgz(x, con, ...)

## S3 method for class 'niftiImage'
io_write_mgz(x, con, ...)

## S3 method for class 'ants.coreANTS_image.ANTsImage'
io_write_mgz(x, con, ...)

## S3 method for class 'array'
io_write_mgz(x, con, vox2ras = NULL, ...)

io_read_nii(
  file,
  method = c("oro", "rnifti", "ants"),
  header_only = FALSE,
  ...
)

io_write_nii(x, con, ...)
```

```
## S3 method for class 'array'
io_write_nii(x, con, vox2ras = NULL, ...)
```

## Arguments

file	file path to read volume data
header_only	whether to read header data only; default is FALSE
format	format of the file to be written; choices are 'auto', 'nifti' or 'mgh'; default is to 'auto' detect the format based on file names, which will save as a 'MGH' file when file extension is 'mgz' or 'mgh', otherwise 'NIfTI' format. We recommend explicitly setting this argument
...	passed to other methods
x	volume data (such as 'NIfTI' image, array, or 'MGH') to be saved
con	file path to store image
vox2ras	a 4x4 transform matrix from voxel indexing (column, row, slice) to scanner (often 'T1-weighted' image) 'RAS' (right-anterior-superior) coordinate
method	method to read the file; choices are 'oro' (using <a href="#">readNIFTI</a> ), 'rnifti' (using <a href="#">readNifti</a> ), and 'ants' (using <a href="#">as_ANTsImage</a> ).

## Format

format of the file; default is auto-detection, other choices are 'nifti' and 'mgh';

## Value

imaging readers return `ieegio_volume` objects.

## Examples

```
library(ieegio)

nifti_file <- "brain.demosubject.nii.gz"

# Use `ieegio_sample_data(nifti_file)`
#   to download sample data


if( ieegio_sample_data(nifti_file, test = TRUE) ) {

  # ---- NIfTI examples -----
  
  file <- ieegio_sample_data(nifti_file)

  # basic read
  vol <- read_volume(file)

  # voxel to scanner RAS
```

```
vol$transforms$vox2ras

# to freesurfer surface
vol$transforms$vox2ras_tkr

# to FSL
vol$transforms$vox2fsl

image(vol$data[,,128], asp = 1, axes = FALSE)

# ---- using other methods -----
# default
vol <- read_volume(file, method = "oro", format = "nifti")
vol$header

# lazy-load nifti
vol2 <- read_volume(file, method = "rnifti", format = "nifti")
vol2$header

# Using ANTsPyx
vol3 <- read_volume(file, method = "ants", format = "nifti")
vol3$header

# ---- write -----
# write as NIfTI
f <- tempfile(fileext = ".nii.gz")

write_volume(vol, f, format = "nifti")

# alternative method
write_volume(vol$header, f, format = "nifti")

# write to mgz/mgh
f2 <- tempfile(fileext = ".mgz")

write_volume(vol, f, format = "mgh")

# clean up
unlink(f)
unlink(f2)

}
```

---

**io\_h5\_valid**

*Check whether a 'HDF5' file can be opened for read/write*

---

**Description**

Check whether a 'HDF5' file can be opened for read/write

**Usage**

```
io_h5_valid(file, mode = c("r", "w"), close_all = FALSE)

io_h5_names(file)
```

**Arguments**

file	path to file
mode	'r' for read access and 'w' for write access
close_all	whether to close all connections or just close current connection; default is false. Set this to TRUE if you want to close all other connections to the file

**Value**

*io\_h5\_valid* returns a logical value indicating whether the file can be opened. *io\_h5\_names* returns a character vector of dataset names.

**Examples**

```
x <- array(1:27, c(3,3,3))
f <- tempfile()

# No data written to the file, hence invalid
io_h5_valid(f, 'r')

io_write_h5(x, f, 'dset')
io_h5_valid(f, 'w')

# Open the file and hold a connection
ptr <- hdf5r::H5File$new(filename = f, mode = 'w')

# Can read, but cannot write
io_h5_valid(f, 'r') # TRUE
io_h5_valid(f, 'w') # FALSE

# However, this can be reset via `close_all=TRUE`
io_h5_valid(f, 'r', close_all = TRUE)
io_h5_valid(f, 'w') # TRUE

# Now the connection is no longer valid
ptr

# clean up
unlink(f)
```

---

**io\_read\_h5***Lazy Load 'HDF5' File via [hdf5r-package](#)*

---

**Description**

Wrapper for class [LazyH5](#), which load data with "lazy" mode - only read part of dataset when needed.

**Usage**

```
io_read_h5(file, name, read_only = TRUE, ram = FALSE, quiet = FALSE)
```

**Arguments**

file	'HDF5' file
name	group/data_name path to dataset (H5D data)
read_only	only used if ram=FALSE, whether the returned <a href="#">LazyH5</a> instance should be read only
ram	load data to memory immediately, default is false
quiet	whether to suppress messages

**Value**

If ram is true, then return data as arrays, otherwise return a [LazyH5](#) instance.

**See Also**

[io\\_write\\_h5](#)

**Examples**

```
file <- tempfile()
x <- array(1:120, dim = c(4,5,6))

# save x to file with name /group/dataset/1
io_write_h5(x, file, '/group/dataset/1', quiet = TRUE)

# read data
y <- io_read_h5(file, '/group/dataset/1', ram = TRUE)
class(y) # array

z <- io_read_h5(file, '/group/dataset/1', ram = FALSE)
class(z) # LazyH5

dim(z)

# clean up
unlink(file)
```

**io\_write\_h5***Save objects to 'HDF5' file without trivial checks***Description**

Save objects to 'HDF5' file without trivial checks

**Usage**

```
io_write_h5(
  x,
  file,
  name,
  chunk = "auto",
  level = 4,
  replace = TRUE,
  new_file = FALSE,
  ctype = NULL,
  quiet = FALSE,
  ...
)
```

**Arguments**

<code>x</code>	an array, a matrix, or a vector
<code>file</code>	path to 'HDF5' file
<code>name</code>	path/name of the data; for example, "group/data_name"
<code>chunk</code>	chunk size
<code>level</code>	compress level from 0 - no compression to 10 - max compression
<code>replace</code>	should data be replaced if exists
<code>new_file</code>	should removing the file if old one exists
<code>ctype</code>	data type such as "character", "integer", or "numeric". If set to NULL then automatically detect types. Note for complex data please store separately the real and imaginary parts.
<code>quiet</code>	whether to suppress messages, default is false
<code>...</code>	passed to other LazyH5\$save

**Value**

Absolute path of the file saved

**See Also**

[io\\_read\\_h5](#)

## Examples

```

file <- tempfile()
x <- array(1:120, dim = 2:5)

# save x to file with name /group/dataset/1
io_write_h5(x, file, '/group/dataset/1', chunk = dim(x))

# load data
y <- io_read_h5(file, '/group/dataset/1')

# read data to memory
y[]

# clean up
unlink(file)

```

LazyH5

*Lazy 'HDF5' file loader*

## Description

provides hybrid data structure for 'HDF5' file

## Public fields

`quiet` whether to suppress messages

## Methods

### Public methods:

- `LazyH5$finalize()`
- `LazyH5$print()`
- `LazyH5$new()`
- `LazyH5$save()`
- `LazyH5$open()`
- `LazyH5$close()`
- `LazyH5$subset()`
- `LazyH5$get_dims()`
- `LazyH5$get_type()`

**Method** `finalize()`: garbage collection method

*Usage:*

`LazyH5$finalize()`

*Returns:* none

**Method** `print()`: overrides print method

*Usage:*

`LazyH5$print()`

*Returns:* self instance

**Method** `new()`: constructor

*Usage:*

`LazyH5$new(file_path, data_name, read_only = FALSE, quiet = FALSE)`

*Arguments:*

`file_path` where data is stored in 'HDF5' format

`data_name` the data stored in the file

`read_only` whether to open the file in read-only mode. It's highly recommended to set this to be true, otherwise the file connection is exclusive.

`quiet` whether to suppress messages, default is false

*Returns:* self instance

**Method** `save()`: save data to a 'HDF5' file

*Usage:*

`LazyH5$save(`

```
x,
chunk = "auto",
level = 7,
replace = TRUE,
new_file = FALSE,
force = TRUE,
ctype = NULL,
size = NULL,
...)
```

*Arguments:*

`x` vector, matrix, or array

`chunk` chunk size, length should matches with data dimension

`level` compress level, from 1 to 9

`replace` if the data exists in the file, replace the file or not

`new_file` remove the whole file if exists before writing?

`force` if you open the file in read-only mode, then saving objects to the file will raise error. Use `force=TRUE` to force write data

`ctype` data type, see [mode](#), usually the data type of `x`. Try `mode(x)` or `storage.mode(x)` as hints.

`size` deprecated, for compatibility issues

`...` passed to self `open()` method

**Method** `open()`: open connection

*Usage:*

```
LazyH5$open(new_dataset = FALSE, robj, ...)
```

*Arguments:*

`new_dataset` only used when the internal pointer is closed, or to write the data  
`robj` data array to save  
... passed to `createDataSet` in `hdf5r` package

**Method** `close()`: close connection

*Usage:*

```
LazyH5$close(all = TRUE)
```

*Arguments:*

`all` whether to close all connections associated to the data file. If true, then all connections, including access from other programs, will be closed

**Method** `subset()`: subset data

*Usage:*

```
LazyH5$subset(..., drop = FALSE, stream = FALSE, envir = parent.frame())
```

*Arguments:*

`drop` whether to apply `drop` the subset  
`stream` whether to read partial data at a time  
`envir` if `i, j, ...` are expressions, where should the expression be evaluated  
`i, j, ...` index along each dimension

*Returns:* subset of data

**Method** `get_dims()`: get data dimension

*Usage:*

```
LazyH5$get_dims(stay_open = TRUE)
```

*Arguments:*

`stay_open` whether to leave the connection opened

*Returns:* dimension of the array

**Method** `get_type()`: get data type

*Usage:*

```
LazyH5$get_type(stay_open = TRUE)
```

*Arguments:*

`stay_open` whether to leave the connection opened

*Returns:* data type, currently only character, integer, raw, double, and complex are available, all other types will yield "unknown"

---

**low-level-read-write** *Low-level file read and write*

---

**Description**

Interfaces to read from or write to files with common formats.

**Usage**

```
io_read_fst(
  con,
  method = c("proxy", "data_table", "data_frame", "header_only"),
  ...,
  old_format = FALSE
)

io_write_fst(x, con, compress = 50, ...)

io_read_ini(con, ...)

io_read_json(con, ...)

io_write_json(
  x,
  con = stdout(),
  ...,
  digits = ceiling(-log10(.Machine$double.eps)),
  pretty = TRUE,
  serialize = TRUE
)

io_read_mat(
  con,
  method = c("auto", "R.matlab", "pymatreader", "mat73"),
  verbose = TRUE,
  on_convert_error = c("warning", "error", "ignore"),
  ...
)

io_write_mat(x, con, method = c("R.matlab", "scipy"), ...)

io_read_yaml(con, ...)

io_write_yaml(x, con, ..., sorted = FALSE)
```

**Arguments**

con	connection or file
-----	--------------------

method	method to read table. For 'fst', the choices are 'proxy' do not read data to memory, query the table when needed; 'data_table' read as <code>data.table</code> ; 'data_frame' read as <code>data.frame</code> ; 'header_only' read 'fst' table header.
	For 'mat', the choices are 'auto' automatically try the native option, and then 'pymatreader' if fails; 'R.matlab' use the native method (provided by <code>readMat</code> ); only support 'MAT 5.0' format; 'pymatreader' use 'Python' library 'pymatreader'; 'mat73' use 'Python' library 'mat73'.
...	passed to internal function calls
old_format	see <code>fst</code>
x	data to write to disk
compress	compress level from 0 to 100; default is 50
digits, pretty	for writing numeric values to 'json' format
serialize	set to TRUE to serialize the data to 'json' format (with the data types, default); or FALSE to save the values without types
verbose	whether to print out the process
on_convert_error	for reading 'mat' files with 'Python' modules, the results will be converted to R objects in the end. Not all objects can be converted. This input defines the behavior when the conversion fails; choices are "error", "warning", or "ignore"
sorted	whether to sort the list; default is FALSE

## Value

The reader functions returns the data extracted from files, mostly as R objects, with few exceptions on some 'Matlab' files. When reading a 'Matlab' file requires using 'Python' modules, `io_read_mat` will try its best effort to convert 'Python' objects to R. However, such conversion might fail. In this case, the result might partially contain 'Python' objects with warnings.

## Examples

```
# ---- fst -----
f <- tempfile(fileext = ".fst")
x <- data.frame(
  a = 1:10,
  b = rnorm(10),
  c = letters[1:10]
)
```

```
io_write_fst(x, con = f)

# default reads in proxy
io_read_fst(f)

# load as data.table
io_read_fst(f, "data_table")

# load as data.frame
io_read_fst(f, "data_frame")

# get header
io_read_fst(f, "header_only")

# clean up
unlink(f)

# ---- json -----
f <- tempfile(fileext = ".json")

x <- list(a = 1L, b = 2.3, c = "a", d = 1+1i)

# default is serialize
io_write_json(x, f)

io_read_json(f)

cat(readLines(f), sep = "\n")

# just values
io_write_json(x, f, serialize = FALSE, pretty = FALSE)

io_read_json(f)

cat(readLines(f), sep = "\n")

# clean up
unlink(f)

# ---- Matlab .mat -----
## Not run:

f <- tempfile(fileext = ".mat")

x <- list(a = 1L, b = 2.3, c = "a", d = 1+1i)

# save as MAT 5.0
io_write_mat(x, f)
```

```

io_read_mat(f)

# require setting up Python environment

io_read_mat(f, method = "pmatreader")

# MAT 7.3 example
sample_data <- ieegio_sample_data("mat_v73.mat")
io_read_mat(sample_data)

# clean up
unlink(f)

## End(Not run)

# ---- yaml -----
f <- tempfile(fileext = ".yaml")

x <- list(a = 1L, b = 2.3, c = "a")
io_write_yaml(x, f)

io_read_yaml(f)

# clean up
unlink(f)

```

**NWBHDF5IO***Creates a NWBHDF5IO file container***Description**

Class definition for 'PyNWB' container; use [read\\_nwb](#) for construction function.

**Active bindings**

`opened` Whether the container is opened.

**Methods****Public methods:**

- [NWBHDF5IO\\$new\(\)](#)
- [NWBHDF5IO\\$get\\_handler\(\)](#)
- [NWBHDF5IO\\$open\(\)](#)

- NWBHDF5IO\$close()
- NWBHDF5IO\$close\_linked\_files()
- NWBHDF5IO\$read()
- NWBHDF5IO\$with()
- NWBHDF5IO\$clone()

**Method new():** Initialize the class

*Usage:*

```
NWBHDF5IO$new(path = NULL, mode = c("r", "w", "r+", "a", "w-", "x"), ...)
```

*Arguments:*

path Path to a '.nwb' file

mode Mode for opening the file

... Other parameters passed to nwb\$NWBHDF5IO

**Method get\_handler():** Get internal file handler. Please make sure you close the handler correctly.

*Usage:*

```
NWBHDF5IO$get_handler()
```

*Returns:* File handler, i.e. 'PyNWB' NWBHDF5IO instance.

**Method open():** Open the connections, must be used together with \$close method. For high-level method, see \$with

*Usage:*

```
NWBHDF5IO$open()
```

*Returns:* container itself

*Examples:*

```
\dontrun{
```

```
# low-level method to open NWB file, for safer methods, see
# `container$with()` below
```

```
container$open()
```

```
data <- container$read()
```

```
# process data...
```

```
# Make sure the container is closed!
container$close()
```

```
}
```

**Method close():** Close the connections (low-level method, see 'with' method below)

*Usage:*

```
NWBHDF5IO$close(close_links = TRUE)
```

*Arguments:*

`close_links` Whether to close all files linked to from this file; default is true

*Returns:* Nothing

**Method** `close_linked_files()`: Close all opened, linked-to files. 'MacOS' and 'Linux' automatically release the linked-to file after the linking file is closed, but 'Windows' does not, which prevents the linked-to file from being deleted or truncated. Use this method to close all opened, linked-to files.

*Usage:*

```
NWBHDF5IO$close_linked_files()
```

*Returns:* Nothing

**Method** `read()`: Read the 'NWB' file from the 'IO' source. Please use along with '\$with' method

*Usage:*

```
NWBHDF5IO$read()
```

*Returns:* 'NWBFile' container

**Method** `with()`: Safe wrapper for reading and handling 'NWB' file. See class examples.

*Usage:*

```
NWBHDF5IO$with(expr, quoted = FALSE, envir = parent.frame())
```

*Arguments:*

`expr` R expression to evaluate

`quoted` Whether `expr` is quoted; default is false

`envir` environment for `expr` to evaluate; default is the parent frame (see `parent.frame()`)

*Returns:* Whatever results generated by `expr`

*Examples:*

```
\dontrun{
```

```
  container$with({
    data <- container$read()
    # process data
  })
}
```

**Method** `clone()`: The objects of this class are cloneable with this method.

*Usage:*

```
NWBHDF5IO$clone(deep = FALSE)
```

*Arguments:*

`deep` Whether to make a deep clone.

## Examples

```

## Not run:

# Running this example requires a .nwb file

library(rnwb)
container <- NWBHDF5IO$new(path = file)
container$with({


  data <- container$read()
  electrode_table <- data$electrodes[convert = TRUE]

})

print(electrode_table)

## End(Not run)

## -----
## Method `NWBHDF5IO$open`
## -----


## Not run:

# low-level method to open NWB file, for safer methods, see
# `container$with()` below

container$open()

data <- container$read()

# process data...

# Make sure the container is closed!
container$close()


## End(Not run)

## -----
## Method `NWBHDF5IO$with`
## -----


## Not run:

container$with({
  data <- container$read()
  # process data
})


```

```
## End(Not run)
```

**plot.ieegio\_surface** *Plot '3D' surface objects*

## Description

Plot '3D' surface objects

## Usage

```
## S3 method for class 'ieegio_surface'
plot(
  x,
  method = c("basic", "full"),
  transform = 1L,
  name = "auto",
  vlim = NULL,
  col = c("black", "white"),
  slice_index = NULL,
  ...
)
```

## Arguments

x	'ieegio_surface' object, see <a href="#">read_surface</a>
method	plot method; 'basic' for just rendering the surfaces; 'full' for rendering with axes and title
transform	which transform to use, can be a 4-by-4 matrix; if the surface contains transform matrix, then this argument can be an integer index of the transform embedded, or the target (transformed) space name; print names(x\$transforms) for choices
name	attribute and name used for colors, options can be 'color' if the surface has color matrix; c('annotations', varname) for rendering colors from annotations with variable varname; c('measurements', varname) for rendering colors from measurements with variable varname; 'time_series' for plotting time series slices; or "flat" for flat color; default is 'auto', which will plot the first available data. More details see 'Examples'.
vlim	when plotting with continuous data (name is measurements or time-series), the value limit used to generate color palette; default is NULL: the range of the values. This argument can be length of 1 ( creating symmetric value range) or 2. If set, then values exceeding the range will be trimmed to the limit
col	color or colors to form the color palette when value data is continuous; when name="flat", the last color will be used
slice_index	when plotting the name="time_series" data, the slice indices to plot; default is to select a maximum of 4 slices
...	ignored

## Examples

```

library(ieegio)

# geometry
geom_file <- "gifti/GzipBase64/sujet01_Lwhite.surf.gii"

# measurements
shape_file <- "gifti/GzipBase64/sujet01_Lwhite.shape.gii"

# time series
ts_file <- "gifti/GzipBase64/fmri_sujet01_Lwhite_projection.time.gii"

if(ieegio_sample_data(geom_file, test = TRUE)) {

  geometry <- read_surface(ieegio_sample_data(geom_file))
  measurement <- read_surface(ieegio_sample_data(shape_file))
  time_series <- read_surface(ieegio_sample_data(ts_file))
  ts_demean <- apply(
    time_series$time_series$value,
    MARGIN = 1L,
    FUN = function(x) {
      x - mean(x)
    }
  )
  time_series$time_series$value <- t(ts_demean)

  # merge measurement & time_series into geometry (optional)
  merged <- merge(geometry, measurement, time_series)
  print(merged)

  # ---- plot method/style -----
  plot(merged, "basic")
  plot(merged, "full")

  # ---- plot data -----
  ## Measurements or annotations

  # the first column of `measurements`
  plot(merged, name = "measurements")

  # equivalent to
  plot(merged, name = list("measurements", 1L))

  # equivalent to
  measurement_names <- names(merged$measurements$data_table)
  plot(merged, name = list("measurements", measurement_names[[1]]))

  ## Time-series
}

```

```
# automatically select 4 slices, trim the color palette
# from -25 to 25
plot(merged, name = "time_series", vlim = c(-25, 25))

plot(
  merged,
  name = "time_series",
  vlim = c(-25, 25),
  slice_index = c(1, 17, 33, 49, 64, 80, 96, 112, 128),
  col = c("#053061", "#2166ac", "#4393c3",
          "#92c5de", "#d1e5f0", "#ffffff",
          "#fddbc7", "#f4a582", "#d6604d",
          "#b2182b", "#67001f")
)
}

}
```

---

plot.ieegio\_volume      *Plot '3D' volume in anatomical slices*

---

## Description

Plot '3D' volume in anatomical slices

## Usage

```
## S3 method for class 'ieegio_volume'
plot(
  x,
  position = c(0, 0, 0),
  center_position = FALSE,
  which = c("coronal", "axial", "sagittal"),
  slice_index = 1L,
  transform = "vox2ras",
  zoom = 1,
  pixel_width = max(zoom/2, 1),
  crosshair_gap = 4,
  crosshair_lty = 2,
  col = c("black", "white"),
  crosshair_col = "#00FF00A0",
  continuous = TRUE,
  vlim = NULL,
  add = FALSE,
  main = "",
```

```

    axes = FALSE,
    background = col[[1]],
    foreground = col[[length(col)]],
    ...,
    .xdata = x$data
)

```

### Arguments

<code>x</code>	'ieegio_volume' object; see <a href="#">read_volume</a>
<code>position</code>	cross-hair focused position
<code>center_position</code>	whether to center canvas at position, default is FALSE
<code>which</code>	which slice to plot; choices are "coronal", "axial", and "sagittal"
<code>slice_index</code>	length of 1: if <code>x</code> has fourth dimension (e.g. 'fMRI'), then which slice index to draw
<code>transform</code>	which transform to apply, can be a 4-by-4 matrix, an integer or name indicating the matrix in <code>x\$transforms</code> ; this needs to be the transform matrix from voxel index to 'RAS' (right-anterior-superior coordinate system), often called 'xform', 'sform', 'qform' in 'NIfTI' terms, or 'Norig' in 'FreeSurfer'
<code>zoom</code>	zoom-in level
<code>pixel_width</code>	pixel size, ranging from 0.05 to 50; default is the half of zoom or 1, whichever is greater; the unit of <code>pixel_width</code> divided by <code>zoom</code> is milliliter
<code>crosshair_gap</code>	the cross-hair gap in milliliter
<code>crosshair_lty</code>	the cross-hair line type
<code>col</code>	color palette for continuous <code>x</code> values
<code>crosshair_col</code>	the cross-hair color; set to NA to hide
<code>continuous</code>	reserved
<code>vlim</code>	the range limit of the data; default is computed from range of <code>x\$data</code> ; data values exceeding the range will be trimmed
<code>add</code>	whether to add the plot to existing underlay; default is FALSE
<code>main, ...</code>	passed to <a href="#">image</a>
<code>axes</code>	whether to draw axes; default is FALSE
<code>background, foreground</code>	background and foreground colors; default is the first and last elements of <code>col</code>
<code>.xdata</code>	default is <code>x\$data</code> , used to speed up the calculation when multiple different angles are to be plotted

### Examples

```

library(ieegio)

nifti_file <- "brain.demosubject.nii.gz"

```

```
# Use `ieegio_sample_data(nifti_file)`  
#   to download sample data  
  
if( ieegio_sample_data(nifti_file, test = TRUE) ) {  
  
  # ---- NIfTI examples -----  
  
  file <- ieegio_sample_data(nifti_file)  
  
  # basic read  
  vol <- read_volume(file)  
  
  par(mfrow = c(1, 3), mar = c(0, 0, 3.1, 0))  
  
  ras_position <- c(-50, -10, 15)  
  
  ras_str <- paste(sprintf("%.0f", ras_position), collapse = ",")  
  
  for(which in c("coronal", "axial", "sagittal")) {  
    plot(x = vol, position = ras_position, crosshair_gap = 10,  
          crosshair_lty = 2, zoom = 3, which = which,  
          main = sprintf("%s T1RAS=[%s]", which, ras_str))  
  }  
  
}  

```

---

pynwb\_module                  *Install 'NWB' via 'pynwb'*

---

## Description

Install 'NWB' via 'pynwb'

## Usage

```
install_pynwb(python_ver = "auto", verbose = TRUE)  
  
pynwb_module(force = FALSE, error_if_missing = TRUE)
```

## Arguments

python_ver	'Python' version, see <a href="#">configure_conda</a> ; default is "auto", which is suggested
verbose	whether to print the installation messages
force	whether to force-reload the module
error_if_missing	whether to raise errors when the module fails to load; default is true

**Value**

A 'Python' module pynwb.

---

read\_bci2000

*Read 'BCI2000' data file*

---

**Description**

Read 'BCI2000' data file

**Usage**

```
read_bci2000(
  file,
  extract_path =getOption("ieegio.extract_path", NULL),
  header_only = FALSE,
  cache_ok = TRUE,
  verbose = TRUE
)
```

**Arguments**

file	file path to the data file
extract_path	location to where the extracted information is to be stored
header_only	whether to only load header data
cache_ok	whether existing cache should be reused; default is TRUE. This input can speed up reading large data files; set to FALSE to delete cache before importing.
verbose	whether to print processing messages; default is TRUE

**Value**

A cached object that is readily to be loaded to memory; see [SignalDataCache](#) for class definition.

**Examples**

```
if( ieegio_sample_data("bci2k.dat", test = TRUE) ) {
  file <- ieegio_sample_data("bci2k.dat")

  x <- read_bci2000(file)
  print(x)

  channel <- x$get_channel(1)

  plot(
    channel$time,
    channel$value,
```

```

    type = "l",
    main = channel$info$Label,
    xlab = "Time",
    ylab = channel$info$Unit
  )
}

```

**read\_brainvis** *Read 'BrainVision' data*

## Description

Read 'BrainVision' data

## Usage

```

read_brainvis(
  file,
  extract_path = getOption("ieegio.extract_path", NULL),
  header_only = FALSE,
  cache_ok = TRUE,
  verbose = TRUE
)

```

## Arguments

<code>file</code>	file path to the data file
<code>extract_path</code>	location to where the extracted information is to be stored
<code>header_only</code>	whether to only load header data
<code>cache_ok</code>	whether existing cache should be reused; default is TRUE. This input can speed up reading large data files; set to FALSE to delete cache before importing.
<code>verbose</code>	whether to print processing messages; default is TRUE

## Value

A cached object that is readily to be loaded to memory; see [SignalDataCache](#) for class definition.

## Examples

```

if( ieegio_sample_data("brainvis.dat", test = TRUE) ) {
  # ensure the header and marker files are downloaded as well
  ieegio_sample_data("brainvis.vhdr")
  ieegio_sample_data("brainvis.dat")
  file <- ieegio_sample_data("brainvis.vmrk")

  x <- read_brainvis(file)
}

```

```

print(x)

x$get_header()

x$get_channel_table()

x$get_annotations()

channel <- x$get_channel(10)

plot(
  channel$time,
  channel$value,
  type = "l",
  main = channel$info$Label,
  xlab = "Time",
  ylab = channel$info$Unit
)
}

```

**read\_edf***Read 'EDF' or 'BDF' data file***Description**

Read 'EDF' or 'BDF' data file

**Usage**

```

read_edf(
  con,
  extract_path = getOption("ieegio.extract_path", NULL),
  header_only = FALSE,
  cache_ok = TRUE,
  begin = 0,
  end = Inf,
  convert = TRUE,
  verbose = TRUE
)

```

**Arguments**

<code>con</code>	file or connection to the data file
<code>extract_path</code>	location to where the extracted information is to be stored
<code>header_only</code>	whether to only load header data
<code>cache_ok</code>	whether existing cache should be reused; default is TRUE. This input can speed up reading large data files; set to FALSE to delete cache before importing.

begin, end	begin and end of the data to read
convert	whether to convert digital numbers to analog signals; default is TRUE
verbose	whether to print processing messages; default is TRUE

**Value**

A cached object that is readily to be loaded to memory; see [SignalDataCache](#) for class definition.

**Examples**

```
# ---- EDF/BDF(+) ----

# Run `ieegio_sample_data("edfPlusD.edf")` to download sample data

# Tun example if the sample data exists
if(ieegio_sample_data("edfPlusD.edf", test = TRUE)) {

  edf_path <- ieegio_sample_data("edfPlusD.edf")

  data <- read_edf(edf_path)

  data$get_header()

  data$get_annotations()

  data$get_channel_table()

  channel <- data$get_channel(1)

  plot(
    channel$time,
    channel$value,
    type = "l",
    main = channel$info$Label,
    xlab = "Time",
    ylab = channel$info$Unit
  )

}
```

**Description**

Read ('BlackRock') 'NEV' 'NSx' data

**Usage**

```
read_nsx(
    file,
    extract_path = getOption("ieegio.extract_path", NULL),
    header_only = FALSE,
    cache_ok = TRUE,
    include_waveform = FALSE,
    verbose = TRUE
)
```

**Arguments**

<code>file</code>	file path to the data file
<code>extract_path</code>	location to where the extracted information is to be stored
<code>header_only</code>	whether to only load header data
<code>cache_ok</code>	whether existing cache should be reused; default is TRUE. This input can speed up reading large data files; set to FALSE to delete cache before importing.
<code>include_waveform</code>	whether to include 'waveform' data (usually for online spike sorting); default is FALSE
<code>verbose</code>	whether to print processing messages; default is TRUE

**Value**

A cached object that is readily to be loaded to memory; see [SignalDataCache](#) for class definition.

`read_nwb`*Read 'NWB' format***Description**

Life cycle: experimental. Read "Neurodata Without Borders" ('NWB' format) file. Unlike other readers `read_nwb` returns low-level 'Python' class handler via `pynwb` module.

**Usage**

```
read_nwb(file, mode = c("r", "w", "r+", "a", "w-", "x"), ...)
```

**Arguments**

<code>file</code>	path to 'NWB' file
<code>mode</code>	file open mode; default is 'r' (read-only)
...	passed to <a href="#">NWBHDF5IO</a> initialize function

**Value**

A NWBHDF5IO instance

**Examples**

```
if(ieegio_sample_data("nwb_sample.nwb", test = TRUE)) {  
    file <- ieegio_sample_data("nwb_sample.nwb")  
  
    # Create NWBIO container  
    container <- read_nwb(file)  
  
    # Open connection  
    container$open()  
  
    # read meta data  
    data <- container$read()  
    data  
  
    # get `test_timeseries` data  
    ts_data <- data$get_acquisition("test_timeseries")  
    ts_data  
  
    # read timeseries data into memory  
    ts_arr <- ts_data$data[]  
    ts_arr  
  
    # Convert Python array to R  
    # using `rpymat::py_to_r(ts_arr)` or  
    as.numeric(ts_arr)  
  
    # Make sure you close the connection  
    container$close()  
  
}  
  
# Requires setting up Python environment  
# run `ieegio::install_pynwb()` to set up environment first  
## Not run:  
  
# Replicating tutorial  
# https://pynwb.readthedocs.io/en/stable/tutorials/general/plot\_file.html  
  
library(rpymat)  
  
# Load Python module  
pynwb <- import("pynwb")  
uuid <- import("uuid")  
datetime <- import("datetime")
```

```

np <- import("numpy")
tz <- import("dateutil.tz")

# 2018L is 2018 as integer
session_start_time <- datetime$datetime(
  2018L, 4L, 25L, 2L, 30L, 3L,
  tzinfo=tz$gettz("US/Pacific"))

# ---- Create NWB file object -----
nwbfile <- pynwb$NWBFfile(
  session_description="Mouse exploring a closed field",
  identifier=py_str(uuid$uuid4()),
  session_start_time=session_start_time,
  session_id="session_4321",
  experimenter=py_list(c("Baggins, Frodo")),
  lab="Bag End Laboratory",
  institution="University of Middle Earth at the Shire",
  experiment_description="Thank you Bilbo Baggins.",
  keywords=py_list(c("behavior", "exploration")))
)

# ---- Add subject -----
subject <- pynwb$file$Subject(
  subject_id="001",
  age="P90D",
  description="mouse 5",
  species="Mus musculus",
  sex="M"
)
nwbfile$subject <- subject

nwbfile

# ---- Add TimeSeries -----
data <- seq(100, 190, by = 10)
time_series_with_rate <- pynwb$TimeSeries(
  name="test_timeseries",
  description="an example time series",
  data=data,
  unit="m",
  starting_time=0.0,
  rate=1.0
)
time_series_with_rate

nwbfile$add_acquisition(time_series_with_rate)

# ---- New Spatial positions -----
position_data <- cbind(
  seq(0, 10, length.out = 50),

```

```

    seq(0, 9, length.out = 50)
  )
position_timestamps = seq(0, 49) / 200

spatial_series_obj = pynwb$behavior$SpatialSeries(
  name="SpatialSeries",
  description="(x,y) position in open field",
  data=position_data,
  timestamps=position_timestamps,
  reference_frame="(0,0) is bottom left corner",
)
spatial_series_obj

position_obj = pynwb$behavior$Position(
  spatial_series=spatial_series_obj)
position_obj

# ---- Behavior Processing Module -----
behavior_module <- nwbfle$create_processing_module(
  name="behavior", description="processed behavioral data"
)
behavior_module$add(position_obj)

nwbfle$processing$behavior

# omit some process

# ---- Write -----
f <- normalizePath(tempfile(fileext = ".nwb"),
  winslash = "/",
  mustWork = FALSE)
io <- pynwb$NWBHDF5IO(f, mode = "w")
io$write(nwbfle)
io$close()

## End(Not run)

```

**Description**

This class is an internal abstract class

**Methods****Public methods:**

- `SignalDataCache$get_header()`
- `SignalDataCache$get_annotations()`
- `SignalDataCache$get_channel_table()`
- `SignalDataCache$get_channel()`
- `SignalDataCache$delete()`

**Method** `get_header():` Get header information, often small list object

*Usage:*

```
SignalDataCache$get_header(...)
```

*Arguments:*

... passed to child methods

**Method** `get_annotations():` Get annotation information, often a large table

*Usage:*

```
SignalDataCache$get_annotations(...)
```

*Arguments:*

... passed to child methods

**Method** `get_channel_table():` Get channel table

*Usage:*

```
SignalDataCache$get_channel_table(...)
```

*Arguments:*

... passed to child methods

**Method** `get_channel():` Get channel data

*Usage:*

```
SignalDataCache$get_channel(x, ...)
```

*Arguments:*

x channel order or label

... passed to child methods

*Returns:* Channel signal with time-stamps inheriting class '`ieegio_get_channel`'

**Method** `delete():` Delete file cache

*Usage:*

```
SignalDataCache$delete(...)
```

*Arguments:*

... passed to child methods

# Index

as\_ANTsImage, 6  
configure\_conda, 25  
data.frame, 15  
data.table, 15  
drop, 13  
fst, 15  
hdf5r-package, 9  
ieegio\_sample\_data, 2  
image, 24  
imaging-surface, 3  
imaging-volume, 5  
install\_pynwb (pynwb\_module), 25  
io\_h5\_names (io\_h5\_valid), 7  
io\_h5\_valid, 7  
io\_read\_fs (imaging-surface), 3  
io\_read\_fst (low-level-read-write), 14  
io\_read\_gii (imaging-surface), 3  
io\_read\_h5, 9, 10  
io\_read\_ini (low-level-read-write), 14  
io\_read\_json (low-level-read-write), 14  
io\_read\_mat (low-level-read-write), 14  
io\_read\_mgz (imaging-volume), 5  
io\_read\_nii (imaging-volume), 5  
io\_read\_yaml (low-level-read-write), 14  
io\_write\_fst (low-level-read-write), 14  
io\_write\_gii (imaging-surface), 3  
io\_write\_h5, 9, 10  
io\_write\_json (low-level-read-write), 14  
io\_write\_mat (low-level-read-write), 14  
io\_write\_mgz (imaging-volume), 5  
io\_write\_nii (imaging-volume), 5  
io\_write\_yaml (low-level-read-write), 14  
  
LazyH5, 9, 11  
low-level-read-write, 14  
mode, 12  
NWBHDF5IO, 17, 30, 31  
plot.ieegio\_surface, 21  
plot.ieegio\_volume, 23  
pynwb\_module, 25  
read.fs.curv, 3  
read.fs.surface, 3  
read\_bci2000, 26  
read\_brainvis, 27  
read\_edf, 28  
read\_nsx, 29  
read\_nwb, 17, 30  
read\_surface, 21  
read\_surface (imaging-surface), 3  
read\_volume, 24  
read\_volume (imaging-volume), 5  
readMat, 15  
readNIfTI, 6  
readNifti, 6  
SignalDataCache, 26, 27, 29, 30, 33  
write\_surface (imaging-surface), 3  
write\_volume (imaging-volume), 5